

Take the pressure off of pouch application




Natura™ Cut-to-Fit Accordion Flange lifts easily to provide generous finger room so coupling is easy and comfortable. Designed for tender abdomens or for anyone who wants everyday comfort.

Natura™ Cut-to-Fit Accordion Flange

NOW Natura™ Cut-to-Fit Accordion Flange

Designed to make pouch application easy and comfortable



Comfortable and discreet:


- In the “UP” position the Accordion Flange provides full finger width space so pouch can be attached with minimal pressure on the abdomen
- In the “DOWN” position it gives the wearer a flat low profile that is extremely flexible and conforms to body contours



“UP” position



“DOWN” position




With the control of a cut-to-fit barrier:

- It's simple to measure and cut to create an opening tailored to each individual patient's stoma size
- Available in Stomahesive® and Duraheasive® adhesive options



Compatible with a range of Natura™ pouch options:

- Pair the Accordion Flange with our state-of-the-art Natura™  comfort pouches in drainable and closed-end options

Ordering Information

Natura™ Accordion Flange: Cut-to-Fit

Product Code	Color Code	Description	Stoma Size	Flange Size	Qty Per Box
421453	●	Duraheasive™ Skin Barrier with Acrylic Collar 1 3/4" (45mm)	1/2" - 3/4" (13-21mm)	1 3/4" (45mm)	5
421454	●	Duraheasive™ Skin Barrier with Acrylic Collar 1 3/4" (45mm)	1/2" - 3/4" (13-21mm)	1 3/4" (45mm)	10
421455	●	Stomahesive™ Skin Barrier with Acrylic Collar 1 3/4" (45mm)	1/2" - 3/4" (13-21mm)	1 3/4" (45mm)	5
421456	●	Stomahesive™ Skin Barrier with Acrylic Collar 1 3/4" (45mm)	1/2" - 3/4" (13-21mm)	1 3/4" (45mm)	10
421457	●	Duraheasive™ Skin Barrier with Acrylic Collar 2 1/4" (57mm)	1/2" - 1 1/4" (13-33mm)	2 1/4" (57mm)	5
421458	●	Duraheasive™ Skin Barrier with Acrylic Collar 2 1/4" (57mm)	1/2" - 1 1/4" (13-33mm)	2 1/4" (57mm)	10
421459	●	Stomahesive™ Skin Barrier with Acrylic Collar 2 1/4" (57mm)	1/2" - 1 1/4" (13-33mm)	2 1/4" (57mm)	5
421460	●	Stomahesive™ Skin Barrier with Acrylic Collar 2 1/4" (57mm)	1/2" - 1 1/4" (13-33mm)	2 1/4" (57mm)	10
421461	●	Duraheasive™ Skin Barrier with Acrylic Collar 2 3/4" (70mm)	1/2" - 1 3/4" (13-45mm)	2 3/4" (70mm)	5
421462	●	Duraheasive™ Skin Barrier with Acrylic Collar 2 3/4" (70mm)	1/2" - 1 3/4" (13-45mm)	2 3/4" (70mm)	10
421463	●	Stomahesive™ Skin Barrier with Acrylic Collar 2 3/4" (70mm)	1/2" - 1 3/4" (13-45mm)	2 3/4" (70mm)	5
421464	●	Stomahesive™ Skin Barrier with Acrylic Collar 2 3/4" (70mm)	1/2" - 1 3/4" (13-45mm)	2 3/4" (70mm)	10

Find out More

1-800-422-8811 Monday – Friday 8:30 am – 7:00 pm

www.ConvaTec.com