


# Purifeed Feeding Tubes



with ENFit™ connector

**Intended use:**

A feeding tube is inserted through the nose or mouth via the oesophagus into the stomach, in order to administer medicine, liquids and/or enteral feeding by means of a syringe or tubing. The product is intended for single use only.



**Product range:**

	Effective length mm	Tube, Outer diameter, mm 			
		1.67	2.0	2.67	3.33
		CH05 REF***	CH06 REF***	CH08 REF***	CH10 REF***
	380	1714503	1714506	1714529	1714531
	480	1714508	1714512	1714515	
<b>X-ray line</b>	980	1714518	1714521	1714524	1714526

\* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.  
 \*\* Languages on labels: English, Swedish, Danish, Norwegian, Finnish.

**Material specifications:**

**Contains latex?** ..... No  
**Contains phthalates?** ..... No  
**Tube** ..... Polyurethane with X-ray line  
**ENFit™ connector.** ..... Polyvinyl chloride  
**ENFit™ connector cap** ..... High density polyethylene  
**ENFit™ connector strap.** ..... Low density polyethylene

**Certifications:**

**CE-marking?** ..... Yes (CE 0086)  
**Product classification** ..... Class IIB (EU)  
**Quality system certified to** ..... ISO 13485  
**Global Medical Device Nomenclature (GMDN)** ..... 14221

**Sterilization:**

**Sterile?** ..... Yes  
**Sterilized by** ..... Ethylene oxide  
**Shelf life** ..... 5 years

**Packaging specifications:**

REF	Primary package	Secondary package/Shelf box		Tertiary package/Shipper			Pallet
	Qty, pcs	Qty, pcs	Dimensions, mm (LxWxD)	Qty, pcs	Dimensions, mm (LxWxD)	Volume, m <sup>3</sup>	Qty, pcs
1714503	1	100	297 x 190 x 366	400	391 x 304 x 743	0.088	6,000
1714506	1	100	297 x 190 x 366	400	391 x 304 x 743	0.088	6,000
1714529	1	100	297 x 190 x 366	400	391 x 304 x 743	0.088	6,000
1714531	1	100	297 x 190 x 366	400	391 x 304 x 743	0.088	6,000
1714508	1	100	297 x 190 x 366	400	391 x 304 x 743	0.088	6,000
1714512	1	100	297 x 190 x 366	400	391 x 304 x 743	0.088	6,000
1714515	1	100	297 x 190 x 366	400	391 x 304 x 743	0.088	6,000
1714518	1	100	297 x 190 x 366	400	391 x 304 x 743	0.088	6,000
1714521	1	100	297 x 190 x 366	400	391 x 304 x 743	0.088	6,000
1714524	1	100	297 x 190 x 366	400	391 x 304 x 743	0.088	6,000
1714526	1	100	297 x 190 x 366	400	391 x 304 x 743	0.088	6,000

# Purifeed Feeding Tubes

with ENFit™ connector



# Unomedical

## EAN codes:

REF	Primary package	Secondary package	Tertiary package
1714503	(01)30768455147237	(01)50768455147231	(01)70768455147235
1714506	(01)30768455147268	(01)50768455147262	(01)70768455147266
1714529	(01)30768455147299	(01)50768455147293	(01)70768455147297
1714531	(01)30768455147312	(01)50768455147316	(01)70768455147310
1714508	(01)30768455147336	(01)50768455147330	(01)70768455147334
1714512	(01)30768455147367	(01)50768455147361	(01)70768455147365
1714515	(01)30768455147381	(01)50768455147385	(01)70768455147389
1714518	(01)30768455147404	(01)50768455147408	(01)70768455147402
1714521	(01)30768455147435	(01)50768455147439	(01)70768455147433
1714524	(01)30768455147466	(01)50768455147460	(01)70768455147464
1714526	(01)30768455147480	(01)50768455147484	(01)70768455147488

## Storage and handling:



**ConvaTec Limited**

First Avenue  
Deeside Industrial Park  
Deeside, Flintshire, CH5 2NU, UK

[www.convatec.com](http://www.convatec.com)



# ConvaTec

For further information on these products, please contact your local ConvaTec sales representative.

With reservations for misprints and alterations in product specifications. Printed version not controlled.

ENFit is a trade mark of GEDSA used with their permission.

All other trade marks and logos are the property of a ConvaTec group company.

© 2018 ConvaTec Inc. AP-019490-MM, ver. OCT 2018 (L1)