

# Intubation Stylets

**Intended use:**

The Intubation Stylet is used as an aid to facilitate insertion of the Endotracheal Tube. The product is intended for single use only.



**List of product range:**

REF*	CH size	Outer diameter, mm	Length, mm
MM64200005	5	1.6	255
MM64200006	6	2.1	275
MM64200010	10	3.3	337
MM64200014	14	4.9	337

**Material specifications:**

**Contains latex?** . . . . . No  
**Contains DEHP?** . . . . . Yes  
**Stylet.** . . . . . Aluminium rod and Polyvinyl chloride tube

**Certification:**

**CE-marking?** . . . . . Yes (CE 0123)  
**Product classification** . . . . . Class IIA (EU)  
**Quality system certified to** . . . . . ISO 13485  
**Global Medical Device Nomenclature (GMDN)** . . . . . 41829

**Sterilisation:**

**Sterile?** . . . . . Yes  
**Sterilized by** . . . . . Ethylene oxide  
**Shelf life** . . . . . 5 years

**Packaging:**

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume, m <sup>3</sup>	
MM64200005	10	505 x 96 x 35	400	413 x 363 x 556	8.4	0.08	400
MM64200006	10	505 x 96 x 35	400	413 x 363 x 556	8.4	0.08	400
MM64200010	10	505 x 96 x 35	400	413 x 363 x 556	8.4	0.08	400
MM64200014	10	505 x 96 x 35	400	413 x 363 x 556	8.4	0.08	400

\* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

**Storage and handling:**



*With reservations for misprints and alterations in product specifications. Printed version not controlled.*

Manufacturer: **Unomedical Sdn Bhd**  
 (15290-P)  
 Address: Bakar Arang Industrial Estate  
 08000 Sungai Petani, Kedah, Malaysia

European Representative: **Unomedical a/s**  
 Address: Kongevejen 2  
 3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.  
 The ConvaTec logo is a trademark of ConvaTec Inc.