


Anesthesia Masks

Intended use:

Anesthesia masks are intended to deliver anesthetic gases to patients. In addition, anesthesia masks are also used to deliver air and/or oxygen during resuscitation and continuous positive airway pressure treatment. The products are intended for single use only.



List of product range:

	REF*	Size	Description
Clear-View Plus Face Masks (with syringe-fill valve)  588UMM 587UMM	589UMM	Adult, Large	Adult, Outer diameter 22 mm, with removable collar ring
	588UMM	Adult, Medium	Adult, Outer diameter 22 mm, with removable collar ring
	587UMM	Adult, Small	Adult, Outer diameter 22 mm, with removable collar ring
	586UMM	Pediatric	Pediatric, Outer diameter 22 mm, with removable collar ring
	585UMM	Infant	Adult, Outer diameter 15 mm
	584UMM	Neonate	Neonate, Outer diameter 15 mm

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
Contains DEHP? Yes
Mask Polyvinyl chloride
Cushion Polyvinyl chloride
Collar ring Polyethylene
Syringe-fill inflation valve Polyvinyl chloride, Polyester, Nitrile, Steel

Certification:

CE-marking? Yes (CE 0120)
Product classification Class IIA (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN). 46232

Sterilisation:

Sterile? No
Shelf life 5 years

Anesthesia Masks

Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume, m ³	
589UMM	*	*	20	381 x 254 x 254	2.99	0.024	500
588UMM	*	*	20	381 x 254 x 279	2.45	0.027	500
587UMM	*	*	20	381 x 254 x 304	1.36	0.029	500
586UMM	*	*	20	381 x 254 x 330	1.63	0.032	1,200
585UMM	*	*	20	279 x 330 x 127	1.09	0.012	500
584UMM	*	*	20	279 x 330 x 152	1.90	0.014	500

* = Single packed - no shelf box

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **Unomedical Inc.**
5701-1S, Ware Road
Address: McAllen, Texas 78503
USA

European Representative: **Unomedical a/s**
Address: Kongevejen 2
3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

