


# Up-Mist™

## Nebulizer & kits

**Intended use:**

Nebulizers are intended for the deposition of particles in order to administer medication via respiration. Up-Mist™ nebulizer is designed for tracheobronchial deposition with a volume of 25 ml. The products are intended for single use only.



**List of product range:**

**Up-Mist™ tracheobrochial nebulizers and kits**

Nebulizer medication capacity: 25 ml

REF*	Item	Size	Description	Oxygen tubing		
				m	Type	Connector
950MM	Up-Mist™ nebulizer		Nebulizer chamber only, no tubing	-	-	-
952MM	Up-Mist™ mask hand held kit	Adult	Nebulizer, T-Adapter, Mouthpiece, with Oxygen tubing	2.1	Standard	Standard
953MM	Up-Mist™ mask kit	Adult	Nebulizer, Aerosol connector, Aerosol mask (elongated), with Oxygen tubing	2.1	Standard	Standard
958MM	Up-Mist™ mask hand held kit	Adult	Nebulizer, T-Adapter, Mouthpiece, no tubing	-	-	-
959MM	Up-Mist™ mask hand held kit	Universal	Nebulizer, T-Adapter, Mouthpiece with 15cm reservoir tube, with Oxygen tubing	2.1	Standard	Standard
3228MM	Up-Mist™ mask kit	Pediatric	Nebulizer, Pediatric Aerosol mask (elongated), with Oxygen tubing	2.1	Standard	Fits-all

\* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

**Material specifications:**

- Contains latex? . . . . . No
- Contains DEHP? . . . . . Yes
- Nebulizer . . . . . Polyethylene
- Aerosol mask . . . . . Polyvinyl chloride
- Mask strap . . . . . Polyvinyl chloride
- Nose clip . . . . . Aluminum
- T-Adapter . . . . . Polypropylene
- Aerosol tubing . . . . . Ethylene vinyl acetate
- Aerosol connector . . . . . High density polyethylene
- Mouthpiece . . . . . High density polyethylene
- Oxygen tubing . . . . . Polyvinyl chloride
- Tube connector . . . . . Polyvinyl chloride

**Certification:**

- CE-marking? . . . . . Yes (CE 0120)
- Product classification . . . . . Class IIA (EU)
- Quality system certified to . . . . . ISO 13485
- Global Medical Device Nomenclature (GMDN) . . . . . 35457

**Sterilisation:**

- Sterile? . . . . . No
- Shelf life . . . . . 5 years

# Up-Mist™

## Nebulizer & kits

### Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m <sup>3</sup>	
950MM	*	*	50	406 x 355 x 254	3.86	0.037	1,250
952MM	*	*	50	355 x 406 x 305	4.76	0.044	1,250
953MM	*	*	50	330 x 482 x 355	7.71	0.056	1,250
958MM	*	*	50	355 x 406 x 305	4.76	0.044	1,250
959MM	*	*	50	381 x 254 x 203	3.63	0.029	1,250
3228MM	*	*	50	355 x 406 x 305	4.76	0.044	1,250

\* = Single packed - no shelf box

\* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

### Storage and handling:



*With reservations for misprints and alterations in product specifications. Printed version not controlled.*

Manufacturer: **Unomedical Inc.**  
5701-1S, Ware Road  
Address: McAllen, Texas 78503  
USA

European Representative: **Unomedical a/s**  
Address: Kongevejen 2  
3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.  
The ConvaTec logo is a trademark of ConvaTec Inc.

