

# A2

## Basic drainage bags

### Intended use:

Urine drainage bag for the collection and monitoring of urine output. The product is intended for single use only.



### List of product range:

#### Drainage bags with non-return valve (NRV)

REF*	Capacity, ml	Tube length, cm
34 001 0-xx	1,500	10
34 003 0-xx	1,500	35
34 005 0-xx	1,500	90
36 030 3-xx	2,000	90
36 105 0-xx	2,000	90

### Material specifications:

**Contains latex?** ..... No  
**Contains DEHP?** ..... Yes  
**Film** ..... Polyvinyl chloride  
**Tube** ..... Polyvinyl chloride  
**Connector** ..... Polyvinyl chloride  
**Cap** . . . Low density polyethylene, High density polyethylene

### Certification:

**CE-marking?** ..... Yes  
**Product classification** ..... Class I (EU)  
**Quality system certified to** ..... ISO 13485  
**Global Medical Device Nomenclature (GMDN)** ..... 58918

### Sterilisation:

**Sterile?** ..... No  
**Shelf life** ..... 5 years

### Packaging:

REF*	Shipping carton				MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume m <sup>3</sup>	
34 001 0-xx	500	375 x 280 x 278	9.0	0.030	24,000
34 003 0-xx	250	375 x 280 x 180	5.8	0.019	9,000
34 005 0-xx	250	388 x 286 x 242	8.0	0.027	14,000
36 030 3-xx	250	375 x 280 x 278	9.8	0.030	12,000
36 105 0-xx	250	388 x 286 x 242	8.7	0.027	14,000

\* Here noted as xx, reference codes contain two additional digits (10-50) representing a language code. Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

### Storage and handling:



*With reservations for misprints and alterations in product specifications. Printed version not controlled.*

Manufacturer: **ConvaTec Limited**  
 First Avenue  
 Address: Deeside Industrial Park  
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.  
The ConvaTec logo is a trademark of ConvaTec Inc.