

# OP-Flex™ : Poole (X-ray visible)

## Suction handles and sets

### Intended use:

Suction handle for intra operative removal of blood and fluids during open surgery. Remove click-on poole for non-dispersed suctioning. Must be attached to a jar/canister via a tube before connecting to a vacuum source of max. 400 mmHg. The product is intended for single use only.



### List of product range:

REF*	Connecting tube Shore A 74		Outer diameter of suction tip mm	Suction handle, length mm	Vacuum control	Number of relieving eyes	CH/FG (handle)		
	Length mm	CH/FG							
OP-Flex™ Poole (X-ray visible)	3407018x	2,100	24	4.6/6.5	175	No	2	20	Medium
	3407118x	3,500	24	4.6/6.5	175	No	2	20	Medium
	3303518x	-	-	4.6/6.5	175	No	2	20	Medium

### Material specifications:

**Contains latex?** . . . . . No  
**Contains DEHP?** . . . . . Yes (340718x, 3407118x)  
**Contains DEHP?** . . . . . No (3303518x)  
**Handle & tip** . . Thermoplastic styrene-butadiene copolymers, Polypropylene  
**Tubing** . . . . . Polyvinyl chloride  
**Connectors** . . . . . Polyvinyl chloride

### Certification:

**CE-marking?** . . . . . Yes (CE 0086)  
**Product classification** . . . . . Class IIA (EU)  
**Quality system certified to** . . . . . ISO 13485  
**Global Medical Device**  
**Nomenclature (GMDN)** . . . . . 35917/36778

### Sterilisation:

**Sterile?** . . . . . Yes  
**Sterilized by** . . . . . Ethylene oxide  
**Shelf life** . . . . . 5 years

### Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m <sup>3</sup>	
3303518x	25	145 x 93 x 446	200	391 x 304 x 455	9.50	0.056	600
3407018x	*	*	50	391 x 304 x 455	7.40	0.056	300
3407118x	*	*	50	391 x 304 x 455	8.50	0.056	150

\* = Single packed - no shelf box

\* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.  
 Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

### Storage and handling:



*With reservations for misprints and alterations in product specifications. Printed version not controlled.*

Manufacturer: **ConvaTec Limited**  
 First Avenue  
 Address: Deeside Industrial Park  
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.  
 The ConvaTec logo is a trademark of ConvaTec Inc.