


# OP-Flex™ : Conventional Yankauer

## Suction handles and sets

**Intended use:**

OP-Flex™ handles and sets are designed for intra operative removal of blood and fluids during open surgery. Must be attached to jar/canister via tube before connecting to vacuum source of max. 400 mmHg. The conventional Yankauers are designed with a standard tip. The product is intended for single use only.



**List of product range:**

	REF*	Connecting tube Shore A 74		Inner/outer diameter of suction, mm	Suction handle, length in mm	Vacuum control	Number of relieving eyes	CH/FG (handle)	
		Length in mm	CH/FG						
OP-Flex™ conventional Yankauer, open tip, single bend	3304718x	-	-	4.6 / 6.5	155	No	4	20	Medium
	3304618x	-	-	4.6 / 6.5	155	Yes	4	20	Medium
	3304818x	-	-	5.0 / 7.5	155	Yes	0	22	Large
	3309318x	-	-	5.0 / 7.5	155	No	4	22	Large
	3409218x	2,100	24	4.6 / 6.5	155	No	4	20	Medium
	3409418x	3,500	24	4.6 / 6.5	155	No	4	20	Medium
	3408918x	2,100	24	4.6 / 6.5	155	Yes	4	20	Medium
	3409118x	3,500	24	4.6 / 6.5	155	Yes	4	20	Medium
	3409518x	4,500	24	4.6 / 6.5	155	No	4	20	Medium
	3409818x	2,100	24	5.0 / 7.5	155	No	0	22	Large
	3409918x	3,500	24	5.0 / 7.5	155	No	0	22	Large
	3409618x	2,100	24	5.0 / 7.5	155	Yes	0	22	Large
	3409718x	3,500	24	5.0 / 7.5	155	Yes	0	22	Large
	3418518x	3,500	24	5.0 / 7.5	155	Yes	4	22	Large
	3418618x	2,100	24	5.0 / 7.5	155	No	4	22	Large
3418718x	3,500	24	5.0 / 7.5	155	No	4	22	Large	

\* Here noted as x, reference codes contain an additional digit (1-5) representing a language code. Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

**Material specifications:**

**Contains latex?** ..... No  
**Contains DEHP** ..... (without connecting tube) No  
 ..... (with connecting tube) Yes  
**Handle & tip** ..... Acrylonitrile-butadiene-styrene resin,  
 Polyvinyl chloride  
**Tubing** ..... Polyvinyl chloride  
**Connectors** ..... Polyvinyl chloride

**Certification:**

**CE-marking?** ..... Yes (CE 0086)  
**Product classification** ..... Class IIA (EU)  
**Quality system certified to** ..... ISO 13485  
**Global Medical Device**  
**Nomenclature (GMDN)** ..... 15917/36778

**Sterilisation:**

**Sterile?** ..... Yes  
**Sterilized by** ..... Ethylene oxide  
**Shelf life** ..... 5 years

# OP-Flex™ : Conventional Yankauer

## Suction handles and sets

### Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m <sup>3</sup>	
3304618x	50	294 x 125 x 360	300	391 x 304 x 743	9.10	0.092	600
3304718x	50	294 x 125 x 360	300	391 x 304 x 743	9.10	0.092	600
3304818x	50	294 x 125 x 360	300	391 x 304 x 743	9.10	0.092	600
3309318x	50	294 x 125 x 360	300	391 x 304 x 743	9.10	0.092	600
3409218x	*	*	50	392 x 306 x 375	6.10	0.045	600
3409418x	*	*	50	392 x 306 x 375	8.50	0.045	150
3408918x	*	*	50	392 x 306 x 375	7.40	0.045	600
3409118x	*	*	50	392 x 306 x 375	8.50	0.045	150
3409518x	*	*	40	392 x 306 x 375	8.50	0.045	120
3409618x	*	*	50	392 x 306 x 375	7.40	0.045	600
3409718x	*	*	50	392 x 306 x 375	8.50	0.045	150
3409818x	*	*	50	392 x 306 x 375	7.40	0.045	600
3409918x	*	*	50	392 x 306 x 375	8.50	0.045	150
3418518x	*	*	50	392 x 306 x 375	8.50	0.045	150
3418618x	*	*	50	392 x 306 x 375	6.10	0.045	600
3418718x	*	*	50	392 x 306 x 375	8.50	0.045	150

\* = Single packed - no shelf box

\* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.

Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

### Storage and handling:



*With reservations for misprints and alterations in product specifications. Printed version not controlled.*

Manufacturer: **ConvaTec Limited**  
First Avenue  
Address: Deeside Industrial Park  
Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.  
The ConvaTec logo is a trademark of ConvaTec Inc.

