


Conventional Yankauer

Suction handles and sets

Intended use:

The Conventional Yankauer is a suction handle intended for the removal of blood and fluids during surgery. The product is intended for single use only.



List of product range:

	REF*	Connecting tube Shore A 74		Inner/outer diameter of suction, mm	Vacuum control	Number of relieving eyes	CH/FG (handle)	
		Length mm	CH/FG					
Conventional Yankauer, open top, single bend (Length of suction handle 210 mm)	3301118x	-	-	4.6/6.5	No	4	20	Medium
	3301018x	-	-	5.0/7.5	No	0	22	Large
	3300618x	-	-	2.5/4.8	No	0	14	Small
	3300718x	-	-	4.1/6.0	No	0	18	Medium
	3400118x	2,100	24	4.6/6.5	No	0	20	Medium
	3400418x	2,100	24	4.6/6.5	No	4	20	Medium
	3400518x	3,000	24	4.6/6.5	No	4	20	Medium
	3400218x	2,100	24	5.0/7.5	No	0	22	Large
	3402718x	2,100	24	5.0/7.5	Yes	0	22	Large
	3400318x	3,000	24	5.0/7.5	No	0	22	Large
	3402818x	3,500	24	5.0/7.5	Yes	0	22	Large
	3401518x	2,100	30	8.0/10.5	No	0	30	X-Large
	3401618x	3,000	30	8.0/10.5	No	0	30	X-Large

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code. Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
 Contains DEHP Yes
 Handle Polyvinyl chloride
 Tubing Polyvinyl chloride
 Other components Polyvinyl chloride

Certification:

CE-marking? Yes (CE 0086)
 Product classification Class IIA (EU)
 Quality system certified to ISO 13485
 Global Medical Device Nomenclature (GMDN) 15917/36778

Sterilisation:

Sterile? Yes
 Sterilized by Ethylene oxide
 Shelf life 5 years

Conventional Yankauer

Suction handles and sets

Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
3300618x	50	294 x 125 x 360	300	391 x 304 x 743	4.60	0.092	600
3300718x	50	294 x 125 x 360	300	391 x 304 x 743	4.55	0.092	600
3301018x	50	294 x 125 x 360	300	391 x 304 x 743	5.60	0.092	600
3301118x	50	294 x 125 x 360	300	391 x 304 x 743	5.60	0.092	600
3400118x	*	*	50	392 x 306 x 340	6.00	0.043	150
3400218x	*	*	50	392 x 306 x 340	6.00	0.043	150
3400318x	*	*	50	392 x 306 x 340	6.00	0.043	150
3400418x	*	*	50	392 x 306 x 340	6.00	0.043	600
3400518x	*	*	50	392 x 306 x 340	6.00	0.043	150
3401518x	*	*	50	392 x 306 x 375	6.70	0.045	150
3401618x	*	*	40	392 x 306 x 375	6.70	0.045	400
3402718x	*	*	50	392 x 306 x 340	6.00	0.043	150
3402818x	*	*	50	392 x 306 x 340	6.00	0.043	150

* = Single packed - no shelf box

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.
Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
First Avenue
Address: Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

