


OP-Flex™ : FilterFlow™

Suction handles and sets

Intended use:

OP-Flex™ FilterFlow™ handles and sets are primarily designed for orthopaedic surgery and combines suctioning function with filtering, removing debris (bone fragments, chips, cement, marrow and blood clots) from the operative area. Must be attached to jar/canister via tube before connecting to a vacuum source.
The product is intended for single use only.



List of product range:

OP-Flex™ FilterFlow™	REF*	Connecting tube Shore A 79		Inner/outer diameter of suction tip mm	Suction tips, length mm	Vacuum control	Number of relieving eyes	CH/FG (handle)	
		Length mm	CH/FG						
Suction set	3418918x	2,700	24	6.0/8.1	150, 230, 280	No	4	24	Large
Suction handle	3309518x	-	-	6.0/8.1	150, 230, 280	No	4	24	Large

Material specifications:

Contains latex? No
Contains DEHP? (Suction set) Yes
 (Suction handle) No
Handle & tip Styrene-acrylonitrile resin
Suction tips Polyvinyl chloride (DEHP-free)
Filter Styrene-acrylonitrile resin
Connecting Tube Polyvinyl chloride
Connectors (Suction set only) Polyvinyl chloride

Certification:

CE-marking? Yes (CE 0086)
Product classification Class IIA (EU)
Quality system certified to ISO 13485
Global Medical Device
Nomenclature (GMDN)
 Suction handle: 35917/Suction set: 36778

Sterilisation:

Sterile? Yes
Sterilized by Ethylene oxide
Shelf life 5 years

Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
3418918x	*	*	30	392 x 306 x 375	6.20	0.045	120
3309518x	20	297 x 190 x 366	40	392 x 306 x 375	4.22	0.045	120

* = Single packed - no shelf box

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.
 Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
 First Avenue
 Address: Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
 The ConvaTec logo is a trademark of ConvaTec Inc.