

Endotracheal Tube

UnoFlex™ Reinforced Endotracheal Tube

Intended use:

Reinforced Endotracheal tubes are intended to be used to assist in the delivery of anaesthetic gases or air to and from the patients' lung. The device is inserted by way of the nose or mouth through the larynx into the trachea for the purpose of maintaining an open airway and mechanically ventilating the patient. It is attached to the anaesthetic machine or respirator via a connector after insertion into the patient. The product is intended for single use only.



List of product range:

UnoFlex™ Reinforced Endotracheal Tubes

Size			Magill, HVLP* Cuff		Magill, HVLP* Cuff - Preloaded with Intubation Stylet		Magill, Plain (without Cuff)		Magill, Plain (without Cuff) - Preloaded with Intubation Stylet	
I.D., mm	O.D., mm	Length, mm	Without Murphy eye	With Murphy eye	Without Murphy eye	With Murphy eye	Without Murphy eye	With Murphy eye	Without Murphy eye	With Murphy eye
			REF**	REF**	REF**	REF**	REF**	REF**	REF**	REF**
2.5	3.9	185					MM61144025	MM61134025	MM61144125	MM61134125
3.0	4.3	190					MM61144030	MM61134030	MM61144130	MM61134130
3.5	4.6	210					MM61144035	MM61134035	MM61144135	MM61134135
4.0	5.6	225					MM61144040	MM61134040	MM61144140	MM61134140
4.5	6.3	255					MM61144045	MM61134045	MM61144145	MM61134145
5.0	7.2	270	MM61224050	MM61214050	MM61224150	MM61214150	MM61144050	MM61134050	MM61144150	MM61134150
5.5	7.7	270	MM61224055	MM61214055	MM61224155	MM61214155	MM61144055	MM61134055	MM61144155	MM61134155
6.0	8.4	295.5	MM61224060	MM61214060	MM61224160	MM61214160	MM61144060	MM61134060	MM61144160	MM61134160
6.5	8.9	315	MM61224065	MM61214065	MM61224165	MM61214165	MM61144065	MM61134065	MM61144165	MM61134165
7.0	9.7	325	MM61224070	MM61214070	MM61224170	MM61214170				
7.5	10.3	335	MM61224075	MM61214075	MM61224175	MM61214175				
8.0	10.8	340	MM61224080	MM61214080	MM61224180	MM61214180				
8.5	11.1	340	MM61224085	MM61214085	MM61224185	MM61214185				
9.0	11.9	340	MM61224090	MM61214090	MM61224190	MM61214190				

*High Volume Low Pressure

** Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
Contains DEHP? Yes
Tube Polyvinyl chloride
Wire Stainless steel
Cuff Polyvinyl chloride
Connector Polypropylene
Pilot balloon Polyvinyl chloride
Inflation line Polyvinyl chloride
Valve Polyvinyl chloride, Nitril, Stainless steel

Certification:

CE-marking? Yes (CE 0123)
Product classification Class IIA(EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN)46967

Sterilisation:

Sterile? Yes
Sterilized by Ethylene oxide
Shelf life 5 years

Endotracheal Tube

UnoFlex™ Reinforced Endotracheal Tube

Packaging:

REF**	Size	Shelf box		Shipping carton				MOQ, pcs.
		Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
MM612240XX*	5.0 to 6.5	10	348 x 150 x 98	100	530 x 325 x 372	5.5	0.06	300
MM612240XX*	7.0 to 9.0	10	390 x 204 x 100	100	530 x 430 x 420	8.2	0.09	300
MM612140XX*	5.0 to 6.0	10	348 x 150 x 98	100	530 x 325 x 372	5.5	0.06	300
MM612140XX*	6.5	10	348 x 150 x 98	100	530 x 325 x 372	5.5	0.06	100
MM612140XX*	7.0 to 9.0	10	390 x 204 x 100	100	530 x 430 x 420	8.2	0.09	300
MM612241XX*	5.0	10	348 x 150 x 98	100	530 x 325 x 372	5.5	0.06	100
MM612241XX*	5.5 to 6.5	10	348 x 150 x 98	100	530 x 325 x 372	5.5	0.06	300
MM612241XX*	7.0 to 9.0	10	390 x 204 x 100	100	530 x 430 x 420	8.2	0.09	300
MM612141XX*	5.0 to 6.5	10	348 x 150 x 98	100	530 x 325 x 372	5.5	0.06	300
MM612141XX*	7.0 to 9.0	10	390 x 204 x 100	100	530 x 430 x 420	8.2	0.09	300
MM611440XX*	2.5 to 6.0	10	348 x 150 x 98	100	530 x 325 x 372	5.5	0.06	300
MM611440XX*	6.5	10	348 x 150 x 98	100	530 x 325 x 372	5.5	0.06	100
MM611340XX*	2.5 to 4.5	10	348 x 150 x 98	100	530 x 325 x 372	5.5	0.06	300
MM611340XX*	5.0	10	348 x 150 x 98	100	530 x 325 x 372	5.5	0.06	100
MM611340XX*	5.5 to 6.5	10	348 x 150 x 98	100	530 x 325 x 372	5.5	0.06	300
MM611441XX*	2.5 to 6.5	10	348 x 150 x 98	100	530 x 325 x 372	5.5	0.06	300
MM611341XX*	2.5 to 6.5	10	348 x 150 x 98	100	530 x 325 x 372	5.5	0.06	100

*XX = size, e.g. 45 = 4.5 mm I.D.

** Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **Unomedical Sdn Bhd**
(15290-P)
Address: Bakar Arang Industrial Estate
08000 Sungai Petani, Kedah, Malaysia

European Representative: **Unomedical a/s**
Address: Kongevejen 2
3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

