

Passive Wound Drainage

Wound drainage system

Intended use:

Passive Wound Drainage is a closed system designed for passive post-operative wound drainage. The passive wound drainage system is connected to a spiral drain which is placed directly into the surgical wound. Spiral drains, with recessed oval eyes in spiral formation, help reducing the risk of tissue contact and growth through opposing eyes. Passive Wound Drainage is primarily used after abdominal surgery in areas with delicate tissue. The product is intended for single use only.



List of product range:

Passive wound drainage system with spiral drain	Drain size (O.D.)	REF*	Tube length, mm	Drain bag capacity, ml
	CH 10	4022518x	900	700
	CH 14	4022618x	900	700
	CH 18	4021818x	900	700
	CH 20	4022918x	900	700
	CH 24	4023018x	900	700
	CH 28	4023118x	900	700

Material specifications:

Contains latex? No
 Contains DEHP? Yes
 Drainage bag Polyvinyl chloride
 Drain and tube Polyvinyl chloride

Certification:

CE-marking? Yes (CE 0086)
 Product classification Class IIA (EU)
 Quality system certified to ISO 13485
 Global Medical Device Nomenclature (GMDN) 35824

Sterilisation:

Sterile? Yes
 Sterilized by Ethylene oxide
 Shelf life 5 years

Packaging:

REF*	Shipping carton				MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume, m ³	
4021818x	50	392 x 306 x 375	1.87	0.045	500
4022518x	50	392 x 306 x 375	1.34	0.045	500
4022618x	50	392 x 306 x 375	1.34	0.045	500
4022918x	50	392 x 306 x 375	4.90	0.045	500
4023018x	50	392 x 306 x 375	5.15	0.045	500
4023118x	50	392 x 306 x 375	5.50	0.045	500

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.
 Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
 First Avenue
 Address: Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
 The ConvaTec logo is a trademark of ConvaTec Inc.