


# FilterFlow™

## Suction handles and sets

### Intended use:

FilterFlow™ suction handles and sets are intended for intra operative removal of blood, fluids and particles, such as bone fragments, chips, cement, marrow and blood clots, during open surgery. The product must be attached to a jar/canister via a tube before connecting to a vacuum source. The products are intended for single use only.



### List of product range:

	REF**	Connecting tube Shore A 78		Inner/outer diameter of suction mm	Length of suction handles* mm	Vacuum control	Number of relieving eyes	CH/FG (handle)	
		Length in mm	CH/FG						
FilterFlow™ suction set	3418818x	2,700	30	5.7/8.1	150, 230, 280	No	4	24	Large
FilterFlow™ suction handle	3309418x	-	-	5.7/8.1	150, 230, 280	No	4	24	Large

\* Delivered with 3 exchangeable suction handles.

### Material specifications:

**Contains latex?** ..... No  
**Contains DEHP** ..... (Set) Yes  
 ..... (Handle) No  
**Handle** ..... Styrene acrylonitrile resin  
**Suction tip** ..... Polyvinyl chloride  
**Filter** ..... Styrene acrylonitrile resin  
**Connecting tube** ..... Polyvinyl chloride  
**Connector** ..... Polyvinyl chloride

### Certification:

**CE-marking?** ..... Yes (CE 0086)  
**Product classification** ..... Class IIA (EU)  
**Quality system certified to** ..... ISO 13485  
**Global Medical Device**  
**Nomenclature (GMDN)** ..... Handle: 35917, Set: 36778

### Sterilisation:

**Sterile?** ..... Yes  
**Sterilized by** ..... Ethylene oxide  
**Shelf life** ..... 5 years

### Packaging:

REF**	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m <sup>3</sup>	
3418818x	*	*	30	392 x 306 x 375	7.31	0.045	120
3309418x	20	297 x 190 x 366	40	392 x 306 x 375	4.22	0.045	120

\* = Single packed - no shelf box

\*\* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.  
 Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

### Storage and handling:



*With reservations for misprints and alterations in product specifications. Printed version not controlled.*

Manufacturer: **ConvaTec Limited**  
 First Avenue  
 Address: Deeside Industrial Park  
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.  
 The ConvaTec logo is a trademark of ConvaTec Inc.