


Siliconized Foley catheter

Female, 2-way, soft valve

Intended use:

The siliconized catheter is a urinary catheter with a silicone fluid treated surface for easy passage through the urethral meatus. The siliconized catheter is intended for short-, medium- and long-term drainage of the bladder. The product is intended for single use only.



List of product range:

Siliconized Foley catheter, Female, 2-way, soft valve

CH size	5-10 ml/cc balloon
12	MM54111210
14	MM54111410
16	MM54111610
18	MM54111810
20	MM54112010
22	MM54112210
24	MM54112410

Material specifications:

Contains latex? Yes
 Contains DEHP? No
 Catheter Latex
 Soft valve Rubber
 Sleeve Rubber
 Catheter surface Silicone emulsion

Certification:

CE-marking? Yes (CE 0123)
 Product classification Class IIB (EU)
 Quality system certified to ISO 13485
 Global Medical Device Nomenclature (GMDN).34917

Sterilisation:

Sterile? Yes
 Sterilized by Gamma radiation
 Shelf life 5 years

Siliconized Foley catheter

Female, 2-way, soft valve

Packaging:

Siliconized Foley catheter, Female, 2-way, soft valve

REF*	Shelf box		Shipping carton				MOQ. pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
MM54111210	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM54111410	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM54111610	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM54111810	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM54112010	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM54112210	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM54112410	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **Unomedical Sdn Bhd**
(15290-P)
Address: Bakar Arang Industrial Estate
08000 Sungai Petani, Kedah, Malaysia

European Representative: **Unomedical a/s**
Address: Kongevejen 2
3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

