


Silicone Drain



Wound drains for Unoversal™ bellows

Intended use:

Silicone drains are placed directly into surgical wounds by means of a trocar or other surgical instrument to perform wound drainage. The product is intended for single use only.



List of product range:

	Drain length mm	Perforation length, mm	Inner/outer diameter of drain, mm 							
			1.6/3.2	2.0/4.0	3.0/5.0	4.0/6.0	4.5/7.0	5.0/8.0	7.0/10.0	8.0/12.0
			CH 10	CH 12	CH 15	CH 18x	CH 21	CH 24	CH 30	CH 36
			REF*	REF*	REF*	REF*	REF*	REF*	REF*	REF*
Silicone drain with spiral perforation with connector for Unoversal™ bellows	500	90	1508718x	1508818x	1508918x	1509018x	1509118x	1509218x	1509318x	1509418x

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code. Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
 Contains DEHP? No
 Drain trocar Silicone
 Drain connector for Unoversal™ Evoprene Rubber
 Coupling Polypropylene

Certification:

CE-marking? Yes (CE 0086)
 Product classification Class IIA (EU)
 Quality system certified to ISO 13485
 Global Medical Device Nomenclature (GMDN). 11305

Sterilisation:

Sterile? Yes
 Sterilized by Ethylene oxide
 Shelf life 5 years

Silicone Drain

Wound drains for Unoversal™ bellows

Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
1508718x	50	360 x 910	200	391 x 304 x 743	7.0	0.092	600
1508818x	50	360 x 910	200	391 x 304 x 743	7.0	0.092	600
1508918x	50	360 x 910	200	391 x 304 x 743	7.5	0.092	600
1509018x	50	360 x 910	150	391 x 304 x 743	7.5	0.092	450
1509118x	50	360 x 910	150	391 x 304 x 743	7.5	0.092	450
1509218x	50	360 x 910	150	391 x 304 x 743	7.5	0.092	450
1509318x	50	360 x 910	150	391 x 304 x 743	7.5	0.092	450
1509418x	50	360 x 910	150	391 x 304 x 743	7.5	0.092	450

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.
Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
First Avenue
Address: Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

