

# Eye-Fix™

## Securement device

### Intended use:

Eye-Fix™ securement device is designed for securing eye pads or protecting anaesthetized eyes against foreign objects. Eye-Fix™ securement device is available with or without a glass frame. The product is intended for single use only.



### List of product range:

	REF*	Size, cm
Eye-Fix™ without glass	640M	10.4 x 9.2
Eye-Fix™ with glass	645M	10.4 x 9.2

### Material specifications:

**Contains latex?** ..... No  
**Contains DEHP** ..... No  
**Carrier** ..... Non-woven polypropylene  
**Adhesive** ..... Medical grade acrylic adhesive  
**Lense (REF Code 645M only)** ..... Styron polystyrene 678E  
**Release liners** ..... 67g/m<sup>2</sup> release paper (one side)

### Certification:

**CE-marking?** ..... Yes  
**Product classification** ..... Class I (EU)  
**Quality system certified to** ..... ISO 13485  
**Global Medical Device Nomenclature (GMDN)** ..... 35715

### Sterilisation:

**Sterile?** ..... No  
**Shelf life** ..... 3 years

### Packaging:

	REF*	Shelf box		Shipping carton				MOQ, pcs.
		Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m <sup>3</sup>	
Eye-Fix™ without glass	640M	100	120 x 110 x 43	1,200	185 x 132 x 370	1.73	0.0088	1,200
Eye-Fix™ with glass	645M	15	120 x 110 x 43	180	185 x 132 x 370	1.03	0.0088	180

\* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

### Storage and handling:



*With reservations for misprints and alterations in product specifications. Printed version not controlled.*

Manufacturer: **ConvaTec Limited**  
 First Avenue  
 Address: Deeside Industrial Park  
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.  
The ConvaTec logo is a trademark of ConvaTec Inc.