

# A3

## Basic drainage bags

**Intended use:**

Urine drainage bag for the collection and monitoring of urine output. The product is intended for single use only.



**List of product range:**

**Drainage bags with non-return valve (NRV)**

| REF*        | Capacity, ml | Tube length, cm | Bottom outlet |
|-------------|--------------|-----------------|---------------|
| 37 205 0-xx | 2,000        | 90              | Push/pull     |
| 37 206 0-xx | 2,000        | 110             | Push/pull     |
| 37 351 0-xx | 2,000        | 90              | Cross valve   |
| 37 452 0-xx | 2,000        | 110             | Cross valve   |

**Material specifications:**

**Contains latex?** . . . . . No  
**Contains DEHP?** . . . . . Yes  
**Bag** . . . . . Polyvinyl chloride  
**Tube** . . . . . Polyvinyl chloride  
**Connector** . . . . . Polyvinyl chloride  
**Cap** . . . Low density polyethylene, High density polyethylene  
**Push/pull valve.** . . Polyvinyl chloride, Low density polyethylene  
**Cross valve.** . . . . . Polypropylene, High density polyethylene

**Certification:**

**CE-marking?** . . . . . Yes  
**Product classification** . . . . . Class I (EU)  
**Quality system certified to** . . . . . ISO 13485  
**Global Medical Device Nomenclature (GMDN)** . . . . . 58922

**Sterilisation:**

**Sterile?** . . . . . No  
**Shelf life** . . . . . 5 years

**Packaging:**

| REF*        | Shipping carton |                 |               |                        | MOQ, pcs. |
|-------------|-----------------|-----------------|---------------|------------------------|-----------|
|             | Qty., pcs.      | Dimensions, mm  | Gross Wt., kg | Volume, m <sup>3</sup> |           |
| 37 205 0-xx | 200             | 388 x 286 x 242 | 8.1           | 0.027                  | 11,200    |
| 37 206 0-xx | 200             | 388 x 286 x 242 | 8.6           | 0.027                  | 9,600     |
| 37 351 0-xx | 150             | 375 x 280 x 278 | 6.2           | 0.029                  | 8,400     |
| 37 452 0-xx | 150             | 375 x 280 x 278 | 6.2           | 0.029                  | 7,200     |

\* Here noted as xx, reference codes contain two additional digits (10-50) representing a language code. Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

**Storage and handling:**



*With reservations for misprints and alterations in product specifications. Printed version not controlled.*

Manufacturer: **ConvaTec Limited**  
 First Avenue  
 Address: Deeside Industrial Park  
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.  
 The ConvaTec logo is a trademark of ConvaTec Inc.