


Suction Catheters













Müly with Fingertip

Intended use:

Suction catheters are intended for oro-, nasopharyngeal and tracheobronchial suctioning. Must be attached to a jar/canister via tube before being connected to vacuum source. The products are intended for single use only.



Product range:

Müly Standard Material (Shore A 78) REF*	Description	CH	O.D. mm	Effective Length cm	Color Coding
Suction catheters, Müly, Straight tip with Fingertip					
0600918x	Müly, Straight, Standard length	5	1.7	465	Grey 
0601118x	Müly, Straight, Standard length	6	2.0	465	Green 
0601318x	Müly, Straight, Standard length	8	2.7	465	Blue 
0601518x	Müly, Straight, Standard length	10	3.3	465	Black 
0601718x	Müly, Straight, Standard length	12	4.0	465	White 
0601918x	Müly, Straight, Standard length	14	4.7	465	Green 
0602118x	Müly, Straight, Standard length	16	5.3	465	Orange 
0602318x	Müly, Straight, Standard length	18	6.0	465	Red 
0601618x	Müly, Straight, Extended length	10	3.3	535	Black 
0601818x	Müly, Straight, Extended length	12	4.0	535	White 
0602018x	Müly, Straight, Extended length	14	4.7	535	Green 
0602218x	Müly, Straight, Extended length	16	5.3	535	Orange 

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code. Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
 Contains DEHP? No
 Suction catheter, fingertip Polyvinyl chloride

Certifications:

CE-marking? Yes (CE 0086)
 Product classification Class IIA (EU)
 Quality system certified to ISO 13485
 Global Medical Device Nomenclature (GMDN) 10749

Sterilization:

Sterile? Yes
 Sterilized by Ethylene oxide
 Shelf life 5 years

Suction Catheters

Mully with Fingertip



Packaging specifications:

REF	Shelf box		Shipping carton				MOQ, pcs
	Qty, pcs	Dimensions, mm (LxWxD)	Qty, pcs	Dimensions, mm (LxWxD)	Gross Wt. kg	Volume m ³	
0600918x	75	145 x 124 x 734	450	391 x 304 x 743	5.7	0.088	900
0601118x	75	145 x 124 x 734	450	391 x 304 x 743	5.9	0.088	900
0601318x	75	145 x 124 x 734	450	391 x 304 x 743	6.3	0.088	900
0601518x	75	145 x 124 x 734	450	391 x 304 x 743	6.8	0.088	900
0601618x	75	145 x 124 x 734	450	391 x 304 x 743	7.2	0.088	900
0601718x	75	145 x 124 x 734	450	391 x 304 x 743	7.2	0.088	900
0601818x	75	145 x 124 x 734	450	391 x 304 x 743	7.7	0.088	900
0601918x	75	145 x 124 x 734	450	391 x 304 x 743	7.7	0.088	900
0602018x	75	145 x 124 x 734	450	391 x 304 x 743	8.1	0.088	900
0602118x	75	145 x 124 x 734	450	391 x 304 x 743	8.1	0.088	900
0602218x	75	145 x 124 x 734	450	391 x 304 x 743	8.6	0.088	900
0602318x	75	145 x 124 x 734	450	391 x 304 x 743	8.5	0.088	900

Storage and handling:



ConvaTec Limited

First Avenue
Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK



ConvaTec Limited
First Avenue, Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK
www.convatec.com

For further information on these products, please contact your local ConvaTec sales representative.

With reservations for misprints and alterations in product specifications. Printed version not controlled.

All trade marks and logos are the property of a ConvaTec group company.

© 2019 ConvaTec Inc. AP-020252-MM, ver. JUN 2019