

A4

Sterile urine bags

Intended use:

The A4 sterile urine bag is intended for the collection of urine. The product is intended for single use only. Fitted with bottom outlet.



List of product range:

Sterile urine bag with non-return valve (NRV), 2000 ml

REF*	Tube length cm	Standard sample port (for use with needle)	Kombikon™ needle-free sample port	Integrated hanger	Bottom outlet / Cross valve
38 103 1S-xx	100	No	no	No	Yes
38 117 1S-xx	110	Yes	no	No	Yes
38 151 1S-xx	120	Yes	No	No	Yes
38 155 1S-xx	120	No	Yes	No	Yes
38 161 1S-xx	110	Yes	No	Yes	Yes
38 162 1S-xx	110	No	Yes	Yes	Yes

* Here noted as xx, reference codes contain two additional digits (10-50) representing a language code. Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

- Contains latex?** No
- Contains DEHP?** Yes
- Film** Polyvinyl chloride
- Tube** Polyvinyl chloride
- Standard sample port** Polyvinyl chloride
- Kombikon™ needle-free sample port** . . . Polyvinyl chloride, Polypropylene, High density polyethylene, polyethylene, Poly(styrene-butadiene-styrene)
- Cap** . . . Low density polyethylene, High density polyethylene
- Hanger** Polypropylene, High density polyethylene

Certification:

- CE-marking?** Yes (CE 0086)
- Product classification** Class I (EU)
- Quality system certified to** ISO 13485
- Global Medical Device Nomenclature (GMDN)**.58921

Sterilisation:

- Sterile?** Yes
- Sterilized by** Ethylene oxide
- Shelf life** 5 years

Packaging:

Sterile urine bag with non-return valve (NRV), 2000 ml

REF*	Shipping carton				MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume, m ³	
38 103 1S-xx	80	375 x 280 x 278	4.8	0.030	3,840
38 117 1S-xx	40	375 x 280 x 278	9.6	0.030	1,920
38 151 1S-xx	50	375 x 280 x 278	5.7	0.030	2,400
38 155 1S-xx	50	375 x 280 x 278	5.7	0.030	2,400
38 161 1S-xx	30	388 x 286 x 242	3.8	0.027	1,440
38 162 1S-xx	30	388 x 286 x 242	3.8	0.027	1,440

* Here noted as xx, reference codes contain two additional digits (10-50) representing a language code.
 Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
 First Avenue
 Address: Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
 The ConvaTec logo is a trademark of ConvaTec Inc.

