


# Silicone Elastomer coated Foley catheter

## Tiemann/Coudé, 2-way, hard valve

**Intended use:**

The Silicone Elastomer coated Foley catheter is a urinary catheter intended for short to medium term drainage of the bladder. For routine and post-operative drainage. The product is intended for single use only.



**List of product range:**

Silicone Elastomer coated Foley catheter, Tiemann/Coudé, 2 way, hard valve

CH size	5-10 ml/cc balloon	30 ml/cc balloon
12	MM51261210	MM51261230
14	MM51261410	MM51261430
16	MM51261610	MM51261630
18	MM51261810	MM51261830
20	MM51262010	MM51262030
22	MM51262210	MM51262230
24	MM51262410	MM51262430
26	MM51262610	MM51262630

**Material specifications:**

Contains latex? . . . . . Yes  
 Contains DEHP? . . . . . Yes  
 Catheter . . . . . Latex  
 Hard valve . . . . . Polypropylene  
 Sleeve . . . . . Polyvinyl chloride  
 Catheter surface . . . . . Silicone Elastomer coating

**Certification:**

CE-marking? . . . . . Yes (CE 0123)  
 Product classification . . . . . Class IIB (EU)  
 Quality system certified to . . . . . ISO 13485  
 Global Medical Device Nomenclature (GMDN). . . . . 34917

**Sterilisation:**

Sterile? . . . . . Yes  
 Sterilized by . . . . . Gamma radiation  
 Shelf life . . . . . 5 years

# Silicone Elastomer coated Foley catheter

Tiemann/Coudé, 2-way, hard valve

## Packaging:

Silicone Elastomer coated Foley catheter, Tiemann/Coudé, 2 way, hard valve

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m <sup>3</sup>	
MM51261210	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51261230	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51261410	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51261430	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51261610	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51261630	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51261810	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51261830	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51262010	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51262030	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51262210	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51262230	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51262410	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51262430	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51262610	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51262630	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800

\* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

## Storage and handling:



*With reservations for misprints and alterations in product specifications. Printed version not controlled.*

Manufacturer: **Unomedical Sdn Bhd**  
(15290-P)  
Address: Bakar Arang Industrial Estate  
08000 Sungai Petani, Kedah, Malaysia

European Representative: **Unomedical a/s**  
Address: Kongevejen 2  
3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.  
The ConvaTec logo is a trademark of ConvaTec Inc.

