


Tracheal Suction Set

Specimen Collectors





Intended use:

Tracheal suction sets are used for the sampling of secretions from the trachea or bronchial tubes for microbiological or cytological examination. For use with suction catheters. The products are intended for single use only.



List of product range:

Trachea Set - Specimen trap

	REF*	Description	
	2400118x	Funnel inlet / male connector outlet, Cap Trachea plug (green)	For use with Vacutip catheters or closed suction systems
	2400618x	Male connector inlet / funnel outlet, Cap Trachea plug (green)	For use with funnel catheters
	2400418x	Male fingertip vacuum control inlet / funnel outlet, Cap Trachea plug (green)	For use with funnel catheters
	2400518x	Male fingertip vacuum control inlet / funnel outlet, Cap, Extra container and cap Trachea plug (green)	For use with funnel catheters

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code. Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
 Contains DEHP? Yes
 Specimen container Polyvinyl chloride
 Funnel inlet Polyvinyl chloride
 Male connector outlet Polyvinyl chloride
 Male connector inlet Polyvinyl chloride
 Male Fingertip vacuum control inlet Polyvinyl chloride
 Funnel outlet Polyvinyl chloride
 Extra container Polystyrene Sterolution
 Extra cap Polypropen
 Identification label Paper

Certification:

CE-marking? Yes (CE 0086)
 Product classification Class I (EU)
 Quality system certified to ISO 13485
 Global Medical Device Nomenclature (GMDN).38472

Sterilisation:

Sterile? Yes
 Sterilized by Ethylene oxide
 Shelf life 5 years

Tracheal Suction Set

Specimen Collectors

Packaging:

REF**	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume, m ³	
2400118x	50	294 x 125 x 360	300	391 x 304 x 743	7.0	0.092	900
2400418x	50	294 x 125 x 360	300	391 x 304 x 743	9.8	0.092	900
2400518x	50	297 x 190 x 366	200	391 x 304 x 743	10.2	0.092	600
2400618x	50	294 x 125 x 360	300	391 x 304 x 743	9.8	0.092	900

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.

Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
First Avenue
Address: Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

