


Spiral drains



Wound drains

Intended use:

A spiral drain is placed directly into surgical wounds by means of a trocar or other surgical instrument to perform wound drainage. Spiral drains, with recessed oval eyes in spiral formation, help reducing the risk of tissue contact and growth through opposing eyes. For removal of blood and body fluids after surgery. Must be connected to a collection device. All drains are radiopaque. The products are intended for single use only.



List of product range:

	Drain length mm	Perforation length, mm	Inner/outer diameter of drain, mm 									
			1.6/2.7	2.0/3.3	2.2/4.0	2.4/4.7	3.0/5.3	3.5/6.0	4.3/6.6	5.3/7.9	6.3/9.2	
			CH 08	CH 10	CH 12	CH 14	CH 16	CH 18	CH 20	CH 24	CH 28	
			REF*	REF*	REF*	REF*	REF*	REF*	REF*	REF*	REF*	REF*
Spiral drain without trocar	500	50	4002118x	4002318x	4002718x	4003118x	4004518x	4004718x	4004918x	4005118x	4005318x	
Spiral drain with solid trocar		60		4006118x	4019218x	4019618x	4019318x	4006518x				

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.
 Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
 Contains DEHP Yes
 Drain trocar Polyvinyl chloride
 Trocar Stainless steel

Certification:

CE-marking? Yes (CE 0086)
 Product classification Class IIA (EU)
 Quality system certified to ISO 13485
 Global Medical Device Nomenclature (GMDN). 11305

Sterilisation:

Sterile? Yes
 Sterilized by Ethylene oxide
 Shelf life 5 years

Spiral drains

Wound drains

Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
4002118x	40	294 x 125 x 360	120	392 x 306 x 375	1.00	0.045	240
4002318x	40	294 x 125 x 360	120	392 x 306 x 375	4.00	0.045	240
4006118x	30	294 x 125 x 360	90	392 x 306 x 375	3.35	0.045	90
4002718x	40	294 x 125 x 360	120	392 x 306 x 375	4.00	0.045	120
4003118x	40	294 x 125 x 360	120	392 x 306 x 375	3.00	0.045	240
4004518x	40	294 x 125 x 360	120	392 x 306 x 375	3.50	0.045	120
4004718x	40	294 x 125 x 360	120	392 x 306 x 375	3.65	0.045	240
4004918x	20	297 x 92 x 239	80	391 x 304 x 243	2.95	0.031	80
4005118x	20	297 x 92 x 239	80	391 x 304 x 243	3.30	0.031	80
4005318x	20	297 x 92 x 239	80	391 x 304 x 243	4.28	0.031	80
4006518x	30	294 x 125 x 360	90	392 x 306 x 375	3.70	0.045	90
4019218x	30	294 x 125 x 360	90	392 x 306 x 375	3.96	0.045	90
4019318x	30	294 x 125 x 360	90	392 x 306 x 375	3.96	0.045	90
4019618x	40	294 x 125 x 360	120	392 x 306 x 375	5.28	0.045	120

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.
Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
First Avenue
Address: Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

