

Oxygen Mask








Medium Concentration

Intended use:

Medium concentration masks are intended for administration of oxygen when precise control of inspired oxygen concentration is not mandatory. During inspiration the patient draws gases both from oxygen flowing into the mask through the tubing, as well as from room air via ports on the sides of the mask. Oxygen concentration delivered is 40%-60% depending on the patient's breathing. The masks are connected directly to a compressed air or oxygen supply. The products are intended for single use only.



Product range:

	REF*	Item	Size	Description	Tube length
	93-103MM 93-103MMV (Vanilla scented mask)	Medium Concentration Oxygen Mask	Adult	Adult Medium Concentration Oxygen Mask (Elongated) with nose clip and elastic strap, No tubing	-
	93-102MM 93-102MMV (Vanilla scented mask)	Medium Concentration Oxygen Mask with tubing	Adult	Adult Medium Concentration Oxygen Mask (Elongated) with nose clip and elastic strap, 2.1m (7') No-crush oxygen tube, Fits-All connector	2.1m (7')
	93-1023MMV (Vanilla scented mask)	Medium Concentration Oxygen Mask with tubing	Adult	Adult Medium Concentration Oxygen Mask (Elongated) with nose clip and elastic strap, 3m (10') No-crush oxygen tube, Fits-All connector	3m (10')
	93-223MM 93-223MMV (Vanilla scented mask)	Medium Concentration Oxygen Mask	Pediatric	Pediatric Medium Concentration Oxygen Mask (Elongated) with nose clip and elastic strap, No tubing	-
	93-222MM 93-222MMV (Vanilla scented mask)	Medium Concentration Oxygen Mask with tubing	Pediatric	Pediatric Medium Concentration Oxygen Mask (Elongated) with nose clip and elastic strap, 2.1m (7') No-crush oxygen tube, Fits-All connector	2.1m (7')
	93-263MM	Medium Concentration Oxygen Mask	Infant	Infant Medium Concentration Oxygen Mask with nose clip and elastic strap, No tubing	-
	93-262MM	Medium Concentration Oxygen Mask with tubing	Infant	Infant Medium Concentration Oxygen Mask with nose clip and elastic strap, 2.1m (7') No-crush oxygen tube, Fits-All connector	2.1m (7')

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
Contains phthalates? No
Contains bisphenols? No
MaskPolyvinyl chloride
Elastic strap Polyester
Nose clip Aluminum
Oxygen nipple connector High density polyethylene
Oxygen tubePolyvinyl chloride
Tube connectorPolyvinyl chloride

Certifications:

CE-marking? Yes (CE0120)
Product classification IIA
Quality system certified to EN ISO 13485:2012

Sterilization:

Sterile? No
Shelf life 5 years

Oxygen Mask

Medium Concentration



Packaging specifications:

REF	Primary package	Secondary package/Polybag	Tertiary package/Shipper			Pallet
	Qty, pcs	Qty, pcs	Qty, pcs	Dimensions, mm (LxWxD)	Gross Wt., kg	Qty, pcs
93-103MM	1	10	50	400 x 300 x 240	2.05	1,600
93-103MMV	1	10	50	400 x 300 x 240	2.05	1,600
93-102MM	1	10	50	400 x 380 x 260	4.00	900
93-102MMV	1	10	50	400 x 380 x 260	4.00	900
93-1023MMV	1	10	50	400 x 380 x 260	4.55	900
93-223MM	1	10	50	350 x 300 x 200	3.47	1,600
93-223MMV	1	10	50	350 x 300 x 200	3.47	1,600
93-222MM	1	10	50	400 x 380 x 260	3.28	1,200
93-222MMV	1	10	50	400 x 380 x 260	3.28	900
93-263MM	1	10	50	350 x 300 x 200	1.20	1,600
93-262MM	1	10	50	400 x 300 x 240	3.47	1,600

EAN codes:

REF	Primary package	Secondary package	Tertiary package
93-103MM	(01)06945099301055	(01)26945099301059	(01)46945099301053
93-103MMV	(01)06945099301307	(01)26945099301301	(01)46945099301305
93-102MM	(01)06945099301062	(01)26945099301066	(01)46945099301060
93-102MMV	(01)06945099301314	(01)26945099301318	(01)46945099301312
93-223MM	(01)06945099301079	(01)26945099301073	(01)46945099301077
93-223MMV	(01)06945099301338	(01)26945099301332	(01)46945099301336
93-1023MMV	(01)06945099301321	(01)26945099301325	(01)46945099301329
93-222MM	(01)06945099301086	(01)26945099301080	(01)46945099301084
93-222MMV	(01)06945099301345	(01)26945099301349	(01)46945099301343
93-263MM	(01)06945099301093	(01)26945099301097	(01)46945099301091
93-262MM	(01)06945099301109	(01)26945099301103	(01)46945099301107

Packaging material:

Primary package	Poly Bag - Low density polyethylene / Insert - 80g Paper
Secondary package	Poly Bag - Low density polyethylene / Sticker - Coated Paper
Tertiary package	Sticker - Coated Paper / Carton - Corrugated Board

Storage and handling:

	FRAGILE - Contents of the transport package are fragile and shall therefore be handled with care
	KEEP AWAY FROM RAIN - Transport package shall be kept away from rain
	THIS WAY UP - Indicates correct upright position of the transport package
	DO NOT USE HAND HOOKS - Hooks are prohibited for handling the transport package



Ningbo Shengyurui Medical Appliances Co., Ltd.

No. 138 Binhaisi Rd.
Hangzhou Bay New Zone, 315336, Ningbo City
Zhejiang Province, P.R. China

Made in China



Shanghai Intl. Holding Corp. GmbH

Eiffestrasse 80, 20537 Hamburg
Germany

Distributed by



ConvaTec Limited
First Avenue, Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

www.convatec.com

For further information on these products, please contact your local ConvaTec sales representative.

With reservations for misprints and alterations in product specifications. Printed version not controlled.

All trade marks and logos are the property of a ConvaTec group company.

© 2017 ConvaTec Inc. AP-017304-MM, ver. FEB 2017