

# Disarp

## Disposable abdominal retracting pad

### Intended use:

The Disarp is an abdominal retraction pad designed to keep soft tissue away from the operative area and absorb blood and/or fluids during surgery. Must be soaked in sterile saline prior to use. The product is intended for single use only.



### List of product range:

Size, mm	REF*
250 x 50 x 50	3404318x
450 x 50 x 50	3404418x

### Material specifications:

Contains latex? ..... No  
 Contains DEHP? ..... No  
 Disarp outside material ..... Polyamide  
 Disarp inside material ..... Polyurethane  
 Disarp nuclear ..... Stainless steel

### Certification:

CE-marking? ..... Yes (CE 0086)  
 Product classification ..... Class IIA (EU)  
 Quality system certified to ..... ISO 13485  
 Global Medical Device Nomenclature (GMDN) ..... 15085

### Sterilisation:

Sterile? ..... Yes  
 Sterilized by ..... Ethylene oxide  
 Shelf life ..... 5 years

### Packaging:

REF*	Shipping carton				MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume, m <sup>3</sup>	
3404318x	25	392 x 306 x 375	1.8	0.045	150
3404418x	25	391 x 304 x 585	2.2	0.072	150

\* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.

Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

### Storage and handling:



*With reservations for misprints and alterations in product specifications. Printed version not controlled.*

Manufacturer: **ConvaTec Limited**  
 First Avenue  
 Address: Deeside Industrial Park  
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.  
 The ConvaTec logo is a trademark of ConvaTec Inc.