


Stomach Tubes

Intended use:

A stomach tube is inserted through the nose or mouth via oesophagus into the stomach. This enables the draining of gastric contents, decompressing the stomach, obtaining a specimen of the gastric contents and administering medicine, liquids and/or enteral feeding. The product is intended for single use only.



List of product range:

CH size	Outer diameter mm	2 eyes		4 eyes	
				Hard introducer/mandrin	Soft introducer/mandrin
		Shore hardness			
		78	78	78	78
		Effective length, mm			
		760	970	970	970
		REF*	REF*	REF*	REF*
10	3.3	2303918x			
12	4.0	2300118x	2304018x		
14	4.7	2300218x		2301418x	2301718x
16	5.3	2300318x		2301518x	2301818x
18	6.0	2300418x		2301618x	2301918x
20	6.7	2302918x	2304118x		
22	7.4	2303018x			
24	8.0	2300718x			
25	8.3	2303118x			
28	9.3	2303218x			
30	10.0	2303318x			
32	10.6	2303418x			
35	11.6	2303518x			

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.
 Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
 Contains DEHP? No
 Tubes Polyvinyl chloride
 Connector Polyvinyl chloride
 Soft introducers Polyethylene
 Hard introducers Polypropylene

Certification:

CE-marking? Yes (CE 0086)
 Product classification Class IIA (EU)
 Quality system certified to ISO 13485
 Global Medical Device Nomenclature (GMDN). 14230

Sterilisation:

Sterile? Yes
 Sterilized by Ethylene oxide
 Shelf life 5 years

Stomach Tubes

Packaging:

REF*	Tube outer diameter, mm	Shelf box		Shipping carton				MOQ, pcs.
		Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
2300118x	4.0	100	297 x 197 x 366	400	391 x 304 x 743	7.2	0.090	800
2300218x	4.7	50	297 x 197 x 366	300	391 x 304 x 743	7.8	0.090	600
2300318x	5.3	50	294 x 125 x 360	300	391 x 304 x 743	7.3	0.090	600
2300418x	6.0	50	294 x 125 x 360	300	391 x 304 x 743	7.9	0.090	600
2300718x	8.0	*	*	100	395 x 309 x 380	8.0	0.046	200
2301418x	4.7	60	297 x 190 x 366	240	391 x 304 x 743	8.0	0.090	480
2301518x	5.3	60	297 x 190 x 366	240	391 x 304 x 743	8.6	0.090	480
2301618x	6.0	60	297 x 190 x 366	240	391 x 304 x 743	9.5	0.090	480
2301718x	4.7	50	294 x 125 x 360	300	391 x 304 x 743	8.0	0.090	600
2301818x	5.3	50	294 x 125 x 360	300	391 x 304 x 743	8.6	0.090	600
2301918x	6.0	50	294 x 125 x 360	300	391 x 304 x 743	9.5	0.090	600
2302918x	6.7	*	*	100	395 x 309 x 380	4.0	0.046	200
2303018x	7.4	*	*	100	395 x 309 x 380	4.1	0.046	200
2303118x	8.3	*	*	100	395 x 309 x 380	4.1	0.046	150
2303218x	9.3	*	*	75	395 x 309 x 380	4.1	0.046	150
2303318x	10.0	*	*	75	395 x 309 x 380	4.6	0.046	150
2303418x	10.6	*	*	75	395 x 309 x 380	4.6	0.046	150
2303518x	11.6	*	*	75	395 x 309 x 380	4.6	0.046	150
2303918x	3.3	100	297 x 190 x 366	400	391 x 304 x 743	7.2	0.090	800
2304018x	4.0	100	297 x 190 x 366	400	391 x 304 x 743	7.2	0.090	800
2304118x	CH20	100	297 x 190 x 366	400	391 x 304 x 743	7.2	0.090	800

* = Single packed - no shelf box

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.
Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
First Avenue
Address: Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

