

# A4 PreKon™

## Closed urine drainage system with siliconized Foley catheter

### Intended use:

Closed urine drainage system with preconnected siliconized Foley catheter for the collection of urine.

The product is intended for single use only.



### List of product range:

#### A4 PreKon™ sterile closed urine drainage system with siliconized Foley catheter

REF*	Catheter size, CH	Item description
43 924 1S-xx	14	Kombikon™ needle-free sample port, 110 cm tubing, non-return valve (NRV), 2000 ml bag capacity, Bottom Outlet/Cross valve, Preconnected to a Siliconized Foley catheter, 2-way, hard valve.
43 926 1S-xx	16	
43 928 1S-xx	18	

\* Here noted as xx, reference codes contain two additional digits (10-50) representing a language code.

Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

### Material specifications:

**Contains latex?** . . . . . Yes  
**Contains DEHP?** . . . . . Yes  
**Tube inlet/outlet, Foil.** . . . . . Polyvinyl chloride  
**KombiKon™ needle-free sample port.** . . . Polyvinyl chloride, Polypropylene, High density polyethylene, Polyethylene, Poly(styrene-butadiene-styrene)  
**Cap** . . . . Low density polyethylene, High density polyethylene  
**Hanger** . . . . . Polypropylene, High density polyethylene  
**Bottom outlet/Cross valve** . . . . . Polypropylene, High density polyethylene  
**Catheter** . . . . . Latex  
**Catheter hard valve** . . . . . Polypropylene  
**Catheter sleeve** . . . . . Polyvinyl chloride  
**Catheter surface** . . . . . Silicone emulsion

### Certification:

**CE-marking?** . . . . . Yes (0086)  
**Product classification** . . . . . Class IIB (EU)  
**Quality system certified to** . . . . . ISO 13485  
**Global Medical Device Nomenclature (GMDN).** . . . . 58921

### Sterilisation:

**Sterile?** . . . . . Yes  
**Sterilized by** . . . . . Ethylene oxide  
**Shelf life** . . . . . 5 years

# A4 PreKon™

## Closed urine drainage system with siliconized Foley catheter

### Packaging:

#### A4 PreKon™ sterile closed urine drainage system with siliconized Foley catheter

REF*	Shipping carton				MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume, m <sup>3</sup>	
43 924 1S-xx	30	582 x 250 x 234	6.95	0.034	2,400
43 926 1S-xx	30	582 x 250 x 234	6.95	0.034	2,400
43 928 1S-xx	30	582 x 250 x 234	6.95	0.034	1,440

\* Here noted as xx, reference codes contain two additional digits (10-50) representing a language code.

Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

### Storage and handling:



*With reservations for misprints and alterations in product specifications. Printed version not controlled.*

Manufacturer: **ConvaTec Limited**  
First Avenue  
Address: Deeside Industrial Park  
Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.  
The ConvaTec logo is a trademark of ConvaTec Inc.

