


OP-Flex™ : Super Flow

Suction handles and sets

Intended use:

OP-Flex™ Super Flow is a suction handle for intra operative removal of blood and fluids during open surgery. Must be attached to a canister via a tube before connecting to a vacuum source of max. 400 mmHg. The product is intended for single use only.



List of product range:

	REF*	Connecting tube Shore A 74		Inner/outer diameter of suction mm	Suction handle, length mm	Vacuum control	Number of relieving eyes	CH/FG (handle)	
		Length mm	CH/FG						
OP-Flex™ Super Flow, open tip, single bend	3405818x	2,020	35	7.5/10.0	180	No	2	30	X-large
	3405918x	3,360	35	7.5/10.0	180	No	2	30	X-large
	3303018x	-	-	7.5/10.0	180	No	2	30	X-large
	3406018x	2,020	35	7.5/10.0	180	Yes	2	30	X-large
	3406218x	3,360	35	7.5/10.0	180	Yes	2	30	X-large
OP-Flex™ Super Flow, open tip, straight	3405518x	3,360	35	7.5/10.0	180	No	2	30	X-large
	3405618x	2,020	35	7.5/10.0	180	Yes	2	30	X-large
	3405718x	3,360	35	7.5/10.0	180	Yes	2	30	X-large
	3302918x	-	-	7.5/10.0	180	Yes	2	30	X-large
Connecting tube for Super Flow	3405218x	2,020	35	-	-	-	-	-	-

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.
Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
Contains DEHP? (with connecting tube) Yes
 (without connecting tube) No
Handle & tip Acrylonitrile-butadiene-styrene resin,
 Polyvinyl chloride
Tubing & connectors Polyvinyl chloride

Certification:

CE-marking? Yes (CE 0086)
Product classification Class IIA (EU)
Quality system certified to ISO 13485
Global Medical Device
Nomenclature (GMDN) 15917/36778

Sterilisation:

Sterile? Yes
Sterilized by Ethylene oxide
Shelf life 5 years

OP-Flex™ : Super Flow

Suction handles and sets

Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
3302918x	50	297 x 190 x 366	100	392 x 306 x 375	4.8	0.045	300
3303018x	50	297 x 190 x 366	100	392 x 306 x 375	4.8	0.045	100
3405218x	*	*	50	392 x 306 x 375	7.9	0.045	150
3405518x	*	*	25	596 x 396 x 205	10.6	0.049	75
3405618x	*	*	25	596 x 396 x 205	9.5	0.049	75
3405718x	*	*	25	596 x 396 x 205	10.5	0.049	75
3405818x	*	*	25	596 x 396 x 205	9.5	0.049	75
3405918x	*	*	25	596 x 396 x 205	10.5	0.049	300
3406018x	*	*	25	596 x 396 x 205	9.5	0.049	300
3406218x	*	*	25	596 x 396 x 205	10.5	0.049	75

* = Single packed - no shelf box

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.
Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
First Avenue
Address: Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

