


Silicone Elastomer coated Foley catheter

Tiemann/Coudé, 2-way, hard valve

Intended use:

The Silicone Elastomer coated Foley catheter is a urinary catheter intended for short to medium term drainage of the bladder. For routine and post-operative drainage. The product is intended for single use only.



List of product range:

Silicone Elastomer coated Foley catheter, Tiemann/Coudé, 2 way, hard valve

CH size	5-10 ml/cc balloon	30 ml/cc balloon
12	MM51261210	MM51261230
14	MM51261410	MM51261430
16	MM51261610	MM51261630
18	MM51261810	MM51261830
20	MM51262010	MM51262030
22	MM51262210	MM51262230
24	MM51262410	MM51262430
26	MM51262610	MM51262630

Material specifications:

Contains latex? Yes
 Contains DEHP? Yes
 Catheter Latex
 Hard valve Polypropylene
 Sleeve Polyvinyl chloride
 Catheter surface Silicone Elastomer coating

Certification:

CE-marking? Yes (CE 0123)
 Product classification Class IIB (EU)
 Quality system certified to ISO 13485
 Global Medical Device Nomenclature (GMDN). 34917

Sterilisation:

Sterile? Yes
 Sterilized by Gamma radiation
 Shelf life 5 years

Silicone Elastomer coated Foley catheter

Tiemann/Coudé, 2-way, hard valve

Packaging:

Silicone Elastomer coated Foley catheter, Tiemann/Coudé, 2 way, hard valve

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
MM51261210	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51261230	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51261410	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51261430	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51261610	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51261630	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51261810	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51261830	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51262010	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51262030	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51262210	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51262230	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51262410	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51262430	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51262610	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800
MM51262630	10	505 x 97 x 35	800	575 x 495 x 522	23	0.16	800

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **Unomedical Sdn Bhd**
(15290-P)
Address: Bakar Arang Industrial Estate
08000 Sungai Petani, Kedah, Malaysia

European Representative: **Unomedical a/s**
Address: Kongevejen 2
3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

