


OP-Flex™: FilterFlow™

Suction handles and sets

Intended use:

OP-Flex™ FilterFlow™ handles and sets are primarily designed for orthopaedic surgery and combines suctioning function with filtering, removing debris (bone fragments, chips, cement, marrow and blood clots) from the operative area. Must be attached to jar/canister via tube before connecting to a vacuum source. The products are intended for single use only.



Product range:

	REF*	Connecting tube Shore A 78			Inner/outer diameter of suction tip mm	Suction tips Length** mm	Vacuum control	Number of relieving eyes	Handle CH/FG
		Length mm	Inner/outer diameter, mm	Funnel CH/FG					
OP-Flex™ FilterFlow™ suction set	3418918x	2,700	5.6/8.0	24	6.0 / 8.1	150, 230, 280	-	4	24
OP-Flex™ FilterFlow™ suction handle	3309518x	-	-	-	6.0 / 8.1	150, 230, 280	-	4	24

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.
 Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.
 ** Delivered with 3 exchangeable suction tips.

Material specifications:

Contains latex?No
Contains DEHP?No
HandleStyrene-acrylonitrile resin
Suction tipPolyvinyl chloride
Filter Styrene-acrylonitrile resin
Connecting tubePolyvinyl chloride
ConnectorPolyvinyl chloride

Certifications:

CE-marking? Yes (CE 0086)
Product classification Class IIA (EU)
Quality system certified to ISO 13485
Global Medical Device Handle: 35917
Nomenclature (GMDN) Set: 36778

Sterilization:

Sterile? Yes
Sterilized by Ethylene oxide
Shelf life 5 years

Packaging specifications:

REF	Shelf box		Shipping carton				MOQ, pcs
	Qty, pcs	Dimensions, mm (LxWxD)	Qty, pcs	Dimensions, mm (LxWxD)	Gross Wt. kg	Volume m³	
3418918x	12	297 x 190 x 366	24	392 x 306 x 375	6.68	0.045	96
3309518x	20	297 x 190 x 366	40	392 x 306 x 375	4.22	0.045	120

Storage and handling:




ConvaTec Limited
 First Avenue
 Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

 **ConvaTec**
 ConvaTec Limited
 First Avenue, Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK
www.convatec.com