


# Redon drains


## Wound drains

### Intended use:

A redon drain is placed directly into surgical wounds by means of a trocar or other surgical instrument to perform wound drainage. A redon drain, with even, circular perforation provides many holes for drainage and permits the surgeon to cut the drain at even angles. Must be connected to a collection device. All drains are radiopaque. The products are intended for single use only.



### List of product range:

	Drain length mm	Perforation length mm	Inner/outer diameter of drain, mm						
			1.6/2.7	2.0/3.3	2.2/4.0	2.4/4.7	3.0/5.3	3.5/6.0	
			CH 08	CH 10	CH 12	CH 14	CH 16	CH 18	
			REF*	REF*	REF*	REF*	REF*	REF*	REF*
Redon drain without trocar	500	70	1500218x	1500418x	1500518x	1500618x	1500718x	1500818x	
		140	1500918x	1501118x	1501218x	1501318x	1501418x	1501518x	

### Material specifications:

Contains latex? ..... No  
 Contains DEHP? ..... Yes  
 Drain tube ..... Polyvinyl chloride

### Certification:

CE-marking? ..... Yes (CE 0086)  
 Product classification ..... Class IIA (EU)  
 Quality system certified to ..... ISO 13485  
 Global Medical Device Nomenclature (GMDN) ..... 11305

### Sterilisation:

Sterile? ..... Yes  
 Sterilized by ..... Ethylene oxide  
 Shelf life ..... 5 years

### Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume m <sup>3</sup>	
1500218x	100	297 x 92 x 239	1,200	391 x 304 x 743	9.1	0.092	1,200
1500418x	100	297 x 92 x 239	1,200	391 x 304 x 743	11.0	0.092	1,200
1500518x	50	297 x 92 x 239	600	391 x 304 x 743	7.9	0.092	600
1500618x	50	297 x 92 x 239	600	391 x 304 x 743	7.9	0.092	600
1500718x	50	297 x 92 x 239	600	391 x 304 x 743	8.5	0.092	600
1500818x	50	297 x 92 x 239	600	391 x 304 x 743	9.1	0.092	600
1500918x	100	297 x 92 x 239	1,200	391 x 304 x 743	9.1	0.092	1,200
1501118x	100	297 x 92 x 239	1,200	391 x 304 x 743	11.0	0.092	1,200
1501218x	50	297 x 92 x 239	600	391 x 304 x 743	7.9	0.092	600
1501318x	50	297 x 92 x 239	600	391 x 304 x 743	7.9	0.092	600
1501418x	50	297 x 92 x 239	600	391 x 304 x 743	8.5	0.092	600
1501518x	50	297 x 92 x 239	600	391 x 304 x 743	9.1	0.092	600

\* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.  
 Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

### Storage and handling:



# Redon drains

Wound drains

*With reservations for misprints and alterations in product specifications. Printed version not controlled.*

Manufacturer: **ConvaTec Limited**  
First Avenue  
Address: Deeside Industrial Park  
Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.  
The ConvaTec logo is a trademark of ConvaTec Inc.

