

UnoMeter™ Safeti™ Plus System

Hourly Diuresis Monitoring System (with Fixed Bag and Luer-Lock Sampling Port)



Intended use:

Urine meter with chamber for accurate measuring of hourly diuresis.
Fitted with bottom outlet. Designed to prevent infection.

List of product range:

REF	
158100150190	2.0L Fixed straps and round hooks, 150 cm, free flow
158100140190	2.0L Fixed straps and round hooks, 110 cm, free flow

Material:

Contains latex **No**
 Contains paper **Yes**
 Contains cotton. **No**
 UnoMeter™ Safeti™ Plus . . . **Polyvinyl Chloride, ABS, Polypropylene, Polyethylene, Styrene Acrylonitrile**
 Airfilter **Teflon, non-woven Polyester backing**
 Luer-Lock. **Polycarbonate, Silicone Rubber**

Standardisation:

CE-marking **CE 0086**
 Quality system certified to. . **ISO 13485**
 GMDN **32072**

Sterilisation:

Sterile **Yes**
 Sterilised by. **Ethylene Oxide**
 Shelf life. **5 years from the date of sterilisation**

UnoMeter™ Safeti™ Plus System

Hourly Diuresis Monitoring System (with Fixed Bag and Luer-Lock Sampling Port)

UnoMeter™ Safeti™ Plus

Packaging:

REF	Shipping carton				MOQ pcs.
	Qty. Pcs.	Dimensions mm (box)	Gross weight kg	Volume M ³	
158100150190	10	400 mm x 359 mm x 291 mm	4.28	0.042	10
158100140190	10	400 mm x 359 mm x 291 mm	4.28	0.042	10

Instruction for use:

Yes.

Storage and handling:

