

Leg bags

Intended use:

Leg bags, intended for collection of urine.
The product is intended for single use only.



List of product range:

Leg bags with non-return valve (NRV)

Straight inlet REF*	Corner inlet REF*	Capacity ml	Tube length cm	Bottom outlet
17 001 0-xx	17 051 0-xx	600	10	Push/pull
17 005 0-xx		600	35	Push/pull
	17 347 0-xx	600	10	Lever tap
17 349 0-xx		600	35	Lever tap
	17 353 0-xx	600	10	Cross valve
28 005 0-xx		750	35	No outlet
29 003 0-xx		750	35	Push/pull
29 005 0-xx		750	10	Push/pull
29 312 0-xx		750	45	Cross valve

* Here noted as xx, reference codes contain two additional digits (10-50) representing a language code.
Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

- Contains latex?** No
- Contains DEHP?** Yes
- Films** Polyvinyl chloride
- Tube** Polyvinyl chloride
- Connector** Polyvinyl chloride
- Cap** Low density polyethylene
- Bottom outlet: Lever tap.** Polyvinyl chloride, Acrylonitrile-butadiene-styrene, Silicone
- Bottom outlet: Push/pull.** Polyvinyl chloride, Low density polyethylene, High density polyethylene
- Bottom outlet: Cross valve.** Polypropylene, High density polyethylene

Certification:

- CE-marking?** Yes
- Product classification** Class I (EU)
- Quality system certified to** ISO 13485
- Global Medical Device Nomenclature (GMDN).**58924

Sterilisation:

- Sterile?** No
- Shelf life** 5 years

Leg bags

Packaging:

Leg bags with non-return valve (NRV)

REF*	Shipping carton				MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume, m ³	
17 001 0-xx	500	375 x 280 x 278	9.55	0.089	6,000
17 005 0-xx	250	375 x 280 x 180	6.15	0.019	18,000
17 347 0-xx	150	375 x 280 x 180	5.04	0.019	10,800
17 349 0-xx	150	375 x 280 x 180	5.54	0.019	7,200
17 051 0-xx	500	375 x 280 x 278	9.55	0.089	28,000
17 353 0-xx	50	375 x 280 x 81	1.40	0.009	6,400
28 005 0-xx	250	375 x 280 x 180	5.45	0.019	18,000
29 003 0-xx	250	382 x 284 x 242	9.90	0.026	14,000
29 005 0-xx	400	382 x 284 x 242	7.80	0.026	22,400
29 312 0-xx	200	375 x 280 x 278	5.40	0.029	9,600

* Here noted as xx, reference codes contain two additional digits (10-50) representing a language code.
Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
First Avenue
Address: Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

