


# Ulmer drains

## Wound drains

**Intended use:**

An Ulmer drain is placed directly into surgical wounds by means of a trocar or other surgical instrument to perform wound drainage. They are intended for the removal of blood and body fluids in wounds after surgery. Ulmer drains have graduated hole sizes for equal suction effect along full drain length. Must be connected to a collection device. All drains are radiopaque. The products are intended for single use only



**List of product range:**

	Drain length mm	Perforation length mm	Inner/outer diameter of drain, mm						
			1.6/2.7	2.0/3.3	2.2/4.0	2.4/4.7	3.0/5.3	3.5/6.0	
			CH 08	CH 10	CH 12	CH 14	CH 16	CH 18	
			REF*	REF*	REF*	REF*	REF*	REF*	REF*
<b>Ulmer drain without trocar</b> 	500	50		1501818x	1501918x				
		60	1502018x						
		80		1502118x					
		100			1502418x				
		112					1502818x	1503018x	
		116					1502618x		
	750	100			1502518x				
		112					1502918x	1503118x	
		116				1502718x			
<b>Ulmer drain with solid trocar</b> 	500	80		1505618x					
		100			1505818x				
		116				1506018x			
	750	80		1505718x					
		100			1505918x				
		112					1506318x	1506518x	
		116				1506118x			
<b>Ulmer drain with connector and clamp for Unoversal drainage system (without trocar)</b> 	500	80		1509518x					
		100			1509618x				
		116				1509718x			

All REF: Shore A hardness: 74

\* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.

Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

**Material specifications:**

Contains latex? . . . . . No  
 Contains DEHP? . . . . . Yes  
 Drain tube. . . . . Polyvinyl chloride  
 Trocar . . . . . Steel

**Certification:**

CE-marking? . . . . . Yes (CE 0086)  
 Product classification . . . . . Class IIA (EU)  
 Quality system certified to . . . . . ISO 13485  
 Global Medical Device Nomenclature (GMDN). . . . . 11305

**Sterilisation:**

Sterile? . . . . . Yes  
 Sterilized by . . . . . Ethylene oxide  
 Shelf life . . . . . 5 years

# Ulmer drains

## Wound drains

### Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m <sup>3</sup>	
1501818x	100	297 x 92 x 239	1,200	391 x 304 x 743	11.0	0.092	1,200
1501918x	50	297 x 92 x 239	600	391 x 304 x 743	7.9	0.092	600
1502018x	100	297 x 92 x 239	1,200	391 x 304 x 743	9.1	0.092	1,200
1502118x	100	297 x 92 x 239	1200	391 x 304 x 743	11.0	0.092	1,200
1502418x	50	297 x 92 x 239	600	391 x 304 x 743	7.9	0.092	600
1502518x	100	297 x 190 x 366	400	391 x 304 x 743	6.5	0.092	400
1502618x	50	297 x 92 x 239	600	391 x 304 x 743	7.9	0.092	600
1502718x	100	297 x 190 x 366	400	391 x 304 x 743	6.5	0.092	400
1502818x	50	297 x 92 x 239	600	391 x 304 x 743	8.5	0.092	600
1502918x	100	297 x 190 x 366	400	391 x 304 x 743	6.5	0.092	400
1503018x	50	297 x 92 x 239	600	391 x 304 x 743	9.1	0.092	600
1503118x	100	297 x 190 x 366	400	391 x 304 x 743	6.5	0.092	400
1505618x	100	297 x 190 x 366	200	392 x 306 x 375	5.8	0.045	200
1505718x	100	297 x 190 x 366	200	392 x 306 x 375	5.8	0.045	200
1505818x	100	297 x 190 x 366	200	392 x 306 x 375	6.4	0.045	200
1505918x	100	297 x 190 x 366	200	392 x 306 x 375	6.4	0.045	200
1506018x	100	297 x 190 x 366	200	392 x 306 x 375	8.0	0.045	200
1506118x	100	297 x 190 x 366	200	392 x 306 x 375	8.0	0.045	200
1506318x	75	297 x 190 x 366	150	392 x 306 x 375	6.4	0.045	150
1506518x	75	297 x 190 x 366	150	392 x 306 x 375	6.9	0.045	150
1509518x	100	297 x 190 x 366	400	391 x 304 x 743	8.1	0.092	800
1509618x	100	297 x 190 x 366	400	391 x 304 x 743	9.0	0.092	1,200
1509718x	100	297 x 190 x 366	400	391 x 304 x 743	10.1	0.092	1,200

\* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.  
 Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

### Storage and handling:



*With reservations for misprints and alterations in product specifications. Printed version not controlled.*

Manufacturer: **ConvaTec Limited**  
 First Avenue  
 Address: Deeside Industrial Park  
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.  
 The ConvaTec logo is a trademark of ConvaTec Inc.

