

Careline™ leg bags

Unsterile

Intended use:

Urine drainage bag with scale for collection of urine output. The device is not intended to have a measurement function. The scale is approximate and for indication purposes only. The product is intended for single use only.



List of product range:

Careline™ leg bags, Unsterile

REF*	Volume ml	Tube length cm	Cut to fit	Bottom outlet	Leg straps	Textile film on backside of bag
12 151 0-xx	500	50	Yes	Cross valve	No	Yes
12 111 8-xx	500	10	No	Cross valve	2 per shelf box	Yes
12 151 8-xx	500	50	Yes	Cross valve	2 per shelf box	Yes
12 511 8-xx	500	10	No	Lever tap	2 per shelf box	Yes
12 551 8-xx	500	50	Yes	Lever tap	2 per shelf box	Yes
24 111 8-xx	750	10	No	Cross valve	2 per shelf box	Yes
24 151 8-xx	750	50	Yes	Cross valve	2 per shelf box	Yes

* Here noted as xx, reference codes contain two additional digits (10-50) representing a language code. Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
Contains DEHP? Yes
Tube inlet/outlet, Foil, Connector Polyvinyl chloride
Cross valve Polypropylene, High density polyethylene
Lever tap Polyvinyl chloride, Acrylonitrile-butadiene-styrene, Silicone
Cap Low density polyethylene
Tube on bottom outlet Silicone
Films Polyvinyl chloride

Certification:

CE-marking? Yes (CE 0086)
Product classification Class I (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN) 58924

Sterilisation:

Sterile? No
Shelf life 5 years

Careline™ leg bags

Unsterile

Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume, m ³	
12 111 8-xx	10	270 x 180 x 56	80	382 x 286 x 242	4.2	0.026	4480
12 151 0-xx	10	270 x 180 x 56	80	382 x 286 x 242	5.8	0.026	8400
12 151 8-xx	10	270 x 180 x 56	80	382 x 284 x 242	5.8	0.026	4480
12 511 8-xx	10	270 x 180 x 56	80	382 x 284 x 242	4.2	0.026	4480
12 551 8-xx	10	270 x 180 x 56	80	382 x 284 x 242	5.8	0.026	4480
24 111 8-xx	10	270 x 180 x 56	80	382 x 284 x 242	4.4	0.026	4480
24 151 8-xx	10	270 x 180 x 56	80	382 x 284 x 242	4.4	0.026	4480

* Here noted as xx, reference codes contain two additional digits (10-50) representing a language code.

Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
First Avenue
Address: Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

