

SE4

Sterile urine bags

Intended use:

Urine drainage bag with scale for collection and monitoring of urine output. Not intended to have a measurement function. Scale is for indication purposes only. Fitted with bottom outlet. The product is intended for single use only



List of product range:

Sterile urine drainage bag with non-return valve (NRV) and sample port

REF*	Capacity, ml	Tube length, cm	Bottom outlet
38 129 1S-xx	2,000	100	Cross valve

Material specifications:

Contains latex? No
Contains DEHP? Yes
Connector with sample port. Polyvinyl chloride, Polypropylene, High density polyethylene/Polyethylene, Poly(styrene-butadiene-styrene)
Cap Low density polyethylene, High density polyethylene
Tube inlet/outlet, foil Polyvinyl chloride
Films Polyvinyl chloride
Cross valve. Polypropylene, High density polyethylene

Certification:

CE-marking? Yes (CE 0086)
Product classification Class I (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN) 58921

Sterilisation:

Sterile? Yes
Sterilized by Ethylene oxide
Shelf life 5 years

Packaging:

REF*	Shipping carton				MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume, m ³	
34 129 1S-xx	65	382 x 378 x 203	5.88	0.029	3,120

* Here noted as xx, reference codes contain additional digits (10-50) representing a language code. Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
 First Avenue
 Address: Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
 The ConvaTec logo is a trademark of ConvaTec Inc.