


Endotracheal Tube

Extra Soft

Intended use:

The Endotracheal Tube is inserted through the mouth or nose into trachea and connected to a ventilation device to secure an airway. The product is intended for single use only.



List of product range:

Endotracheal tube, Extra Soft, Magill shape without cuff

Without Murphy Eye REF*	With Murphy Eye REF*	Size	I.D., mm	O.D., mm	Length, mm	ISO 15 mm connector
MM63140020	MM63130020	2.0	2.0	2.9	150	Clear
MM63140025	MM63130025	2.5	2.5	3.4	175	Clear
MM63140030	MM63130030	3.0	3.0	4.1	185	Red
MM63140035	MM63130035	3.5	3.5	4.8	185	Red
MM63140040	MM63130040	4.0	4.0	5.4	205	Green
MM63140045	MM63130045	4.5	4.5	6.1	225	Green
MM63140050	MM63130050	5.0	5.0	6.7	245	Blue
MM63140055	MM63130055	5.5	5.5	7.3	275	Blue
MM63140060	MM63130060	6.0	6.0	8.0	285	Blue

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Product specific information:

- Colour coded connector
- Standard graduation marks
- Murphy eyes
- Standard Magill curve

Material specifications:

- Contains latex? No
- Contains DEHP? Yes
- Connector Polypropylene
- Tube Polyvinyl chloride

Certification:

- CE-marking? Yes (CE 0123)
- Product classification Class IIA (EU)
- Quality system certified to ISO 13485
- Global Medical Device Nomenclature (GMDN).46967

Sterilisation:

- Sterile? Yes
- Sterilized by Ethylene oxide
- Shelf life 5 years

Endotracheal Tube

Extra Soft

Packaging:

Endotracheal tube, Extra Soft, Magill shape without cuff

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume, m ³	
MM63140XXX	10	348 x 150 x 98	100	530 x 325 x 372	6.2	0.07	300
MM63130XXX	10	348 x 150 x 98	100	530 x 325 x 372	6.2	0.07	300**

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

**Size 2.5 MOQ=100

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **Unomedical Sdn Bhd**
(15290-P)
Address: Bakar Arang Industrial Estate
08000 Sungai Petani, Kedah, Malaysia

European Representative: **Unomedical a/s**
Address: Kongevejen 2
3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

