


UnoFlex™ Reinforced Tracheal Tubes















Intended use:

Tracheal Tubes are inserted through the mouth or nose into the trachea and connected to ventilation equipment to secure an airway. The products are intended for single use only.











Product range:

UnoFlex™ Reinforced Tracheal Tubes with HVLP* cuff

With Murphy eye REF	Without Murphy eye REF	With stylet and Murphy eye REF	Size		Length mm*	CRD mm**	Connector colour***
			I.D. mm	O.D. mm			
UM61214030	UM61224030	UM61214130	3.0	5.3	175	12	Red 
UM61214035	UM61224035	UM61214135	3.5	5.8	195	12	Red 
UM61214040	UM61224040	UM61214140	4.0	6.3	215	14	Green 
UM61214045	UM61224045	UM61214145	4.5	6.8	240	14	Green 
UM61214050	UM61224050	UM61214150	5.0	7.3	256	17	Blue 
UM61214055	UM61224055	UM61214155	5.5	7.8	286	17	Blue 
UM61214060	UM61224060	UM61214160	6.0	8.4	295	20	Blue 
UM61214065	UM61224065	UM61214165	6.5	9.0	301	20	Black 
UM61214070	UM61224070	UM61214170	7.0	9.6	315	25	Black 
UM61214075	UM61224075	UM61214175	7.5	10.2	325	25	White 
UM61214080	UM61224080	UM61214180	8.0	10.8	337	26	White 
UM61214085	UM61224085	UM61214185	8.5	11.3	337	26	Green 
UM61214090	UM61224090	UM61214190	9.0	11.9	337	28	Green 
UM61214095	UM61224095	UM61214195	9.5	12.4	337	28	Orange 

UnoFlex™ Reinforced Tracheal Tubes without HVLP* cuff

With Murphy eye REF	Without Murphy eye REF	With stylet and Murphy eye REF	Size		Length mm*	Connector colour***
			I.D. mm	O.D. mm		
UM61134030	UM61144030	UM61134130	3.0	5.3	175	Red 
UM61134035	UM61144035	UM61134135	3.5	5.8	195	Red 
UM61134040	UM61144040	UM61134140	4.0	6.3	215	Green 
UM61134045	UM61144045	UM61134145	4.5	6.8	240	Green 
UM61134050	UM61144050	UM61134150	5.0	7.3	256	Blue 
UM61134055	UM61144055	UM61134155	5.5	7.8	286	Blue 
UM61134060	UM61144060	UM61134160	6.0	8.4	295	Blue 
UM61134065	UM61144065	UM61134165	6.5	9.0	301	Black 

* HVLP = High Volume Low Pressure; Total length of tube excluding connector
 ** CRD = Cuff Resting Diameter
 *** 15 mm ISO connector (ISO 8836)

Material specifications:

Contains latex?	No
Contains phthalates?	No
Tube	Polyvinyl chloride
Wire	Stainless steel
Cuff	Polyvinyl chloride
Inflation line	Polyvinyl chloride
X-ray line	Polyvinyl chloride
Connector	Polypropylene
Pilot balloon	Polyvinyl chloride
Valve	Polyvinyl chloride, Polypropylene, Acrylonitrile-butadiene-styrene, Silicone, Stainless steel
Stylet	Aluminum
Sheath	Polyvinyl chloride

Certifications:

CE-marking?	Yes (0123)
Product classification	Class IIA (EU)
Quality system certified to	ISO 13485
Global Medical Device Nomenclature (GMDN)46569

Sterilization:

Sterile?	Yes
Sterilized by	Ethylene oxide
Shelf life	5 years

Packaging specifications:

REF*	Sizes	Primary package/ Single peelpack		Secondary package/ Shelf box		Tertiary package/ Shipper				Pallet Qty pcs
		Qty pcs	Dimensions, mm (LxW)	Qty pcs	Dimensions, mm (LxWxD)	Qty pcs	Dimensions, mm (LxWxD)	Gross Wt. kg	Volume m ³	
UnoFlex™ Reinforced Tracheal Tubes with HVLP cuff										
UM612140XX UM612240XX	3.0 to 4.5	1	274 x 91	10	280 x 125 x 90	100	470 x 295 x 270	5.0	0.037	1,800
UM612140XX UM612240XX	5.0 to 9.5	1	347 x 91	10	350 x 150 x 95	100	495 x 370 x 320	6.0	0.059	1,200
UnoFlex™ Reinforced Tracheal Tubes without HVLP cuff										
UM611340XX UM611440XX	3.0 to 4.5	1	274 x 91	10	280 x 125 x 90	100	470 x 295 x 270	5.0	0.037	1,800
UM611340XX UM611440XX	5.0 to 6.5	1	347 x 91	10	350 x 150 x 95	100	495 x 370 x 320	6.0	0.059	1,200
UnoFlex™ Reinforced Tracheal Tubes with preloaded stylet										
UM612141XX	3.0 to 4.0	1	274 x 91	10	280 x 125 x 90	100	470 x 295 x 270	5.0	0.037	1,800
UM611341XX	3.0 to 4.0	1	274 x 91	10	280 x 125 x 90	100	470 x 295 x 270	5.0	0.037	1,800
UM612141XX	4.5 to 9.5	1	347 x 91	10	350 x 150 x 95	100	495 x 370 x 320	6.0	0.059	1,200
UM611341XX	4.5 to 6.5	1	347 x 91	10	350 x 150 x 95	100	495 x 370 x 320	6.0	0.059	1,200

*XX = size, e.g. 45 = 4.5 mm I.D.

EAN codes:

UnoFlex™ Reinforced Tracheal Tubes with HVLP cuff and Murphy eye

REF	Primary package	Secondary package	Tertiary package
UM61214030	(01)30768455151425	(01)50768455151429	(01)70768455151423
UM61214035	(01)30768455151432	(01)50768455151436	(01)70768455151430
UM61214040	(01)30768455151449	(01)50768455151443	(01)70768455151447
UM61214045	(01)30768455151456	(01)50768455151450	(01)70768455151454
UM61214050	(01)30768455151463	(01)50768455151467	(01)70768455151461
UM61214055	(01)30768455151470	(01)50768455151474	(01)70768455151478
UM61214060	(01)30768455151487	(01)50768455151481	(01)70768455151485
UM61214065	(01)30768455151494	(01)50768455151498	(01)70768455151492
UM61214070	(01)30768455151500	(01)50768455151504	(01)70768455151508
UM61214075	(01)30768455151517	(01)50768455151511	(01)70768455151515
UM61214080	(01)30768455151524	(01)50768455151528	(01)70768455151522
UM61214085	(01)30768455151531	(01)50768455151535	(01)70768455151539
UM61214090	(01)30768455151548	(01)50768455151542	(01)70768455151546
UM61214095	(01)30768455151555	(01)50768455151559	(01)70768455151553

UnoFlex™ Reinforced Tracheal Tubes with HVLP cuff, without Murphy eye

REF	Primary package	Secondary package	Tertiary package
UM61224030	(01)30768455151289	(01)50768455151283	(01)70768455151287
UM61224035	(01)30768455151296	(01)50768455151290	(01)70768455151294
UM61224040	(01)30768455151302	(01)50768455151306	(01)70768455151300
UM61224045	(01)30768455151319	(01)50768455151313	(01)70768455151317
UM61224050	(01)30768455151326	(01)50768455151320	(01)70768455151324
UM61224055	(01)30768455151333	(01)50768455151337	(01)70768455151331
UM61224060	(01)30768455151340	(01)50768455151344	(01)70768455151348
UM61224065	(01)30768455151357	(01)50768455151351	(01)70768455151355
UM61224070	(01)30768455151364	(01)50768455151368	(01)70768455151362
UM61224075	(01)30768455151371	(01)50768455151375	(01)70768455151379
UM61224080	(01)30768455151388	(01)50768455151382	(01)70768455151386
UM61224085	(01)30768455151395	(01)50768455151399	(01)70768455151393
UM61224090	(01)30768455151401	(01)50768455151405	(01)70768455151409
UM61224095	(01)30768455151418	(01)50768455151412	(01)70768455151416

UnoFlex™ Reinforced Tracheal Tubes without HVLP cuff, with Murphy eye

REF	Primary package	Secondary package	Tertiary package
UM61134030	(01)30768455154198	(01)50768455154192	(01)70768455154196
UM61134035	(01)30768455154204	(01)50768455154208	(01)70768455154202
UM61134040	(01)30768455154211	(01)50768455154215	(01)70768455154219
UM61134045	(01)30768455154228	(01)50768455154222	(01)70768455154226
UM61134050	(01)30768455154235	(01)50768455154239	(01)70768455154233
UM61134055	(01)30768455154242	(01)50768455154246	(01)70768455154240
UM61134060	(01)30768455154259	(01)50768455154253	(01)70768455154257
UM61134065	(01)30768455154266	(01)50768455154260	(01)70768455154264

EAN codes:

UnoFlex™ Reinforced Tracheal Tubes without HVLP cuff and Murphy eye

REF	Primary package	Secondary package	Tertiary package
UM61144030	(01)30768455154099	(01)50768455154093	(01)70768455154097
UM61144035	(01)30768455154105	(01)50768455154109	(01)70768455154103
UM61144040	(01)30768455154112	(01)50768455154116	(01)70768455154110
UM61144045	(01)30768455154129	(01)50768455154123	(01)70768455154127
UM61144050	(01)30768455154136	(01)50768455154130	(01)70768455154134
UM61144055	(01)30768455154143	(01)50768455154147	(01)70768455154141
UM61144060	(01)30768455154150	(01)50768455154154	(01)70768455154158
UM61144065	(01)30768455154167	(01)50768455154161	(01)70768455154165

UnoFlex™ Reinforced Tracheal Tubes with preloaded stylet, HVLP cuff and Murphy eye

REF	Primary package	Secondary package	Tertiary package
UM61214130	(01)30768455154426	(01)50768455154420	(01)70768455154424
UM61214135	(01)30768455154433	(01)50768455154437	(01)70768455154431
UM61214140	(01)30768455154440	(01)50768455154444	(01)70768455154448
UM61214145	(01)30768455154457	(01)50768455154451	(01)70768455154455
UM61214150	(01)30768455154464	(01)50768455154468	(01)70768455154462
UM61214155	(01)30768455154471	(01)50768455154475	(01)70768455154479
UM61214160	(01)30768455154488	(01)50768455154482	(01)70768455154486
UM61214165	(01)30768455154495	(01)50768455154499	(01)70768455154493
UM61214170	(01)30768455154501	(01)50768455154505	(01)70768455154509
UM61214175	(01)30768455154518	(01)50768455154512	(01)70768455154516
UM61214180	(01)30768455154525	(01)50768455154529	(01)70768455154523
UM61214185	(01)30768455154532	(01)50768455154536	(01)70768455154530
UM61214190	(01)30768455154549	(01)50768455154543	(01)70768455154547
UM61214195	(01)30768455154556	(01)50768455154550	(01)70768455154554





UnoFlex™ Reinforced Tracheal Tubes with preloaded stylet and Murphy eye, without HVLP cuff

REF	Primary package	Secondary package	Tertiary package
UM61134130	(01)30768455154587	(01)50768455154581	(01)70768455154585
UM61134135	(01)30768455154594	(01)50768455154598	(01)70768455154592
UM61134140	(01)30768455154600	(01)50768455154604	(01)70768455154608
UM61134145	(01)30768455154617	(01)50768455154611	(01)70768455154615
UM61134150	(01)30768455154624	(01)50768455154628	(01)70768455154622
UM61134155	(01)30768455154631	(01)50768455154635	(01)70768455154639
UM61134160	(01)30768455154648	(01)50768455154642	(01)70768455154646
UM61134165	(01)30768455154655	(01)50768455154659	(01)70768455154653

Packaging material:

Primary package	Peelback pouch - Uncoated paper, XPP film
Secondary package	Ivory cardboard
Tertiary package	Corrugated cardboard

Storage and handling:

	FRAGILE - Contents of the transport package are fragile and shall therefore be handled with care
	KEEP DRY
	THIS WAY UP - Indicates correct upright position of the transport package
	DO NOT USE HAND HOOKS - Hooks are prohibited for handling the transport package



Well Lead Medical Co.,Ltd.
 C-4 Jinhui Industrial Estate
 Hualong, Panyu, Guangzhou
 511434 P.R.China

Made in China



Shanghai Intl. Holding Corp. GmbH
 Eiffestr. 80, 20537 Hamburg
 Germany

Distributed by



ConvaTec Limited
 First Avenue, Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK
www.convatec.com

For further information on these products, please contact your local ConvaTec sales representative.

With reservations for misprints and alterations in product specifications. Printed version not controlled.
 All trade marks and logos are the property of a ConvaTec group company.

© 2020 ConvaTec Inc. AP-019166-MM, ver. NOV 2020