


Frazier

Suction handle

Intended use:

Frazier suction handle, for removal of blood during surgery requiring extremely precise suctioning. The product has a stylet to keep the tip open whilst bent in shape. The product is intended for single use only.



List of product range:

| | CH 06 REF | CH 08 REF | CH 10 REF | CH 12 REF |
|---|--------------|--------------|--------------|--------------|
| Frazier suction handle with long tip (OP-Flex™ system)* | | 309.02.008L | 309.02.010L | 309.02.012L |
| Frazier suction handle* | 309.02.006 | 309.02.008 | 309.02.010 | 309.02.012 |

* Picture: Top - Frazier suction handle with long tip, Bottom - Frazier suction handle.

Material specifications:

Contains latex? No
Contains DEHP? Yes
Tip Alloy of aluminum
Transparent handle Polymethyl methacrylate
White handle Acrylonitrile-butadiene-styrene
Stylet Polyvinyl chloride with X-ray line

Certification:

CE-marking? Yes (CE 0123)
Product classification Class IIA (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN) 35917

Sterilisation:

Sterile? Yes
Sterilized by Ethylene oxide
Shelf life 5 years

Packaging:

| REF** | Shipping carton | | | | MOQ. pcs. |
|-------------|-----------------|----------------|---------------|------------------------|-----------|
| | Qty., pcs. | Dimensions, mm | Gross Wt., kg | Volume, m ³ | |
| 309.02.006 | 50 | 310 x 470 x 90 | 0.556 | 0.013 | 25 |
| 309.02.008 | 50 | 310 x 470 x 90 | 0.618 | 0.013 | 50 |
| 309.02.008L | 25 | 310 x 470 x 90 | 0.669 | 0.013 | 25 |
| 309.02.010 | 50 | 310 x 470 x 90 | 0.669 | 0.013 | 50 |
| 309.02.010L | 25 | 310 x 470 x 90 | 0.665 | 0.013 | 25 |
| 309.02.012 | 50 | 310 x 470 x 90 | 0.665 | 0.013 | 50 |
| 309.02.012L | 25 | 310 x 470 x 90 | 0.665 | 0.013 | 25 |

**Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **Unomedical Sdn Bhd**
 (15290-P)
 Address: Bakar Arang Industrial Estate
 08000 Sungai Petani, Kedah, Malaysia

European Representative: **Unomedical a/s**
 Address: Kongevejen 2
 3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
 The ConvaTec logo is a trademark of ConvaTec Inc.