


Cath-Fix™

Securement device

Intended use:

Cath-Fix™ securement device is designed for securement of catheters, tubes and cables. The product is intended for single use only.



List of product range:

	REF*	Size, cm
Cath-Fix™, Small	630M-S	2.6 x 3.65
Cath-Fix™, Medium	630M-M	5.2 x 6.50
Cath-Fix™, Large	630M-L	8.1 x 10.3

Material specifications:

Contains latex? No
Contains DEHP? No
Carrier Non-woven polypropylene
Adhesive Medical grade acrylic
Release liner Release paper 67g/m² (one side)

Certification:

CE-marking? Yes
Product classification Class I (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN) 34864

Sterilisation:

Sterile? No
Shelf life 3 years

Packaging:

	REF*	Shelf box		Shipping carton				MOQ, pcs.
		Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
Cath-Fix™, Small	630M-S	100	93 x 68 x 39	1,600	420 x 140 x 100	0.77	0.006	1,600
Cath-Fix™, Medium	630M-M	100	93 x 68 x 39	1,600	426 x 156 x 110	1.27	0.006	1,600
Cath-Fix™, Large	630M-L	100	106 x 106 x 41	1,200	185 x 132 x 370	1.70	0.009	1,200

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: ConvaTec Limited
 First Avenue
Address: Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.