

Niko-Fix™

Securement device

Intended use:

Niko-Fix™ securement device is a sterile device for short term securement of peripheral intravenous lines.

The product is intended for single use only.



List of product range:

	REF*	Size, cm
Niko-Fix™	NF-92M	7.2 x 8.5

Material specifications:

Contains latex? No
Contains DEHP? No
Carrier Non-woven, spun laced material
Adhesive Medical grade acrylic
Cushion pad Thermal bonded viscose-Rayon Polypropylene
Anchoring of pad Perforated double adhesive fastape
Backing paper One side silicone coated kraft
Pouch (paper) Grid lacquered medical kraft
Pouch (film) Polyester and polythene

Certification:

CE-marking? Yes (0086)
Product classification Class I (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN) 56631

Sterilisation:

Sterile? Yes
Sterilized by Gamma radiation
Shelf life 3 years

Packaging:

	REF*	Shelf box		Shipping carton				MOQ, pcs.
		Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m³	
Niko-Fix™	NF-92M	100	121 x 105 x 152	1,000	560 x 155 x 258	4.14	0.0224	40,000

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
 First Avenue
 Address: Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
 The ConvaTec logo is a trademark of ConvaTec Inc.