

Suction Catheters

Müly with fingertip

Intended use:

Suction catheters are intended for oro-, nasopharyngeal and tracheobronchial suctioning. Must be attached to a jar/canister via tube before being connected to vacuum source.
The product is intended for single use only.



List of product range:

	Müly Standard Material (Shore A 78) REF**	Item description	CH	O.D. mm	Length* mm	Color coding
Suction catheters, Müly, straight tip, with fingertip						
	0600918x	Müly, straight, Standard length	5	1.7	465	Grey 
	0601118x	Müly, straight, Standard length	6	2.0	465	Green 
	0601318x	Müly, straight, Standard length	8	2.3	465	Blue 
	0601518x	Müly, straight, Standard length	10	2.7	465	Black 
	0601718x	Müly, straight, Standard length	12	3.3	465	White 
	0601918x	Müly, straight, Standard length	14	4.0	465	Green 
	0602118x	Müly, straight, Standard length	16	4.7	465	Orange 
	0602318x	Müly, straight, Standard length	18	5.3	465	Red 
	0601618x	Müly, straight, Extended length	10	6.0	535	Black 
	0601818x	Müly, straight, Extended length	12	2.7	535	White 
	0602018x	Müly, straight, Extended length	14	3.3	535	Green 
	0602218x	Müly, straight, Extended length	16	4.0	535	Orange 

* Effective length

** Here noted as x, reference codes contain an additional digit (1-5) representing a language code.

Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
Contains DEHP? Yes
Suction catheter, fingertip Polyvinyl chloride

Certification:

CE-marking? Yes (CE 0086)
Product classification Class IIA (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN) 10749

Sterilisation:

Sterile? Yes
Sterilized by Ethylene oxide
Shelf life 5 years

Suction Catheters

Mülly with fingertip

Packaging:

Suction catheters, Mülly, straight tip, with fingertip

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume, m³	
Suction catheters, Müllý, straight tip, with fingertip							
0600918x	100	145 x 124 x 734	600	391 x 304 x 743	8.0	0.088	1,200
0601118x	100	145 x 124 x 734	600	391 x 304 x 743	8.0	0.088	1,200
0601318x	100	145 x 124 x 734	600	391 x 304 x 743	8.7	0.088	1,200
0601518x	100	145 x 124 x 734	600	391 x 304 x 743	8.7	0.088	1,200
0601618x	100	145 x 124 x 734	600	391 x 304 x 743	8.8	0.088	1,200
0601718x	100	145 x 124 x 734	600	391 x 304 x 743	8.9	0.088	1,200
0601818x	100	145 x 124 x 734	600	391 x 304 x 743	9.4	0.088	1,200
0601918x	100	145 x 124 x 734	600	391 x 304 x 743	9.0	0.088	1,200
0602018x	100	145 x 124 x 734	600	391 x 304 x 743	10.5	0.088	1,200
0602118x	100	145 x 124 x 734	600	391 x 304 x 743	9.8	0.088	1,200
0602218x	100	145 x 124 x 734	600	391 x 304 x 743	11.2	0.088	1,200
0602318x	100	145 x 124 x 734	600	391 x 304 x 743	10.1	0.088	1,200

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.

Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **ConvaTec Limited**
First Avenue
Address: Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

