

# Aerosol Masks, Tracheostomy Masks, Face Tent for mask nebulization


**Intended use:**

Aerosol masks, tracheostomy masks and face tents are intended for nebulization via a mask or face tent in order to administer medication via respiration. The product is intended for single use only.




**List of product range:**


**Aerosol masks, Elongated (under the chin) with adjustable nose-clip and elastic braid\***

	REF**	Size	Description	Length of aerosol tubing, m
100MM	100MM	Adult	Aerosol mask (elongated) only	-
	118MM	Adult	Aerosol mask (elongated), with aerosol tubing	1.8
	220MM	Pediatric	Aerosol mask (elongated) only	-
220MM	221MM	Pediatric	Aerosol mask (elongated), with Aerosol tubing	1.5


**Aerosol masks, Standard (under the lip) with adjustable nose-clip and elastic braid\***

	REF**	Size	Description	Length of aerosol tubing, m
260MM	200MM	Adult	Aerosol mask (Standard) only	-
	260MM	Infant	Aerosol mask (Standard) only	-

**Trach masks with adjustable elastic braid\***

	REF**	Size	Description	Length of aerosol tubing, m
245MM	245MM	Adult	Trach mask with swivel connector	-
246MM	246MM	Pediatric	Trach mask with swivel connector	-

**Face tents**

	REF**	Size	Description	Length of aerosol tubing, m
110MM	110MM	Adult	Face tent	-

\* Connectors accept standard 22 mm aerosol tubing

\*\* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

**Material specifications:**

Contains latex? . . . . . No  
 Contains DEHP? . . . . . Yes  
 Aerosol mask . . . . . Polyvinyl chloride  
 Trach mask . . . . . Polyvinyl chloride  
 Face tent . . . . . Polyvinyl chloride  
 Nose clip . . . . . Aluminum  
 Elastic braid . . . . . Polyester (Latex-Free)  
 Aerosol tubing . . . . . Ethylene vinyl acetate  
 Aerosol connector . . . . . High density polyethylene

**Certification:**

CE-marking? . . . . . Yes (CE 0120)  
 Product classification . . . . . Class IIA (EU)  
 Quality system certified to . . . . . ISO 13485  
 Global Medical Device  
 Nomenclature (GMDN) . . . . . 35142/35179

**Sterilisation:**

Sterile? . . . . . No  
 Shelf life . . . . . 5 years

# Aerosol Masks, Tracheostomy Masks, Face Tent for mask nebulization

## Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m <sup>3</sup>	
100MM	*	*	50	381 x 254 x 203	2.04	0.019	1,250
110MM	*	*	50	406 x 355 x 254	4.54	0.037	1,250
118MM	*	*	50	406 x 609 x 584	7.26	0.144	1,250
200MM	*	*	50	381 x 254 x 203	1.81	0.019	1,250
220MM	*	*	50	381 x 254 x 203	1.81	0.019	1,250
221MM	*	*	50	406 x 609 x 584	4.99	1.440	1,250
245MM	*	*	50	381 x 254 x 203	2.04	0.019	1,250
246MM	*	*	50	381 x 254 x 203	2.04	0.019	1,250
260MM	*	*	50	381 x 254 x 203	1.36	0.019	1,250

\* = Single packed - no shelf box

\* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

## Storage and handling:



*With reservations for misprints and alterations in product specifications. Printed version not controlled.*

Manufacturer: **Unomedical Inc.**  
5701-1S, Ware Road  
Address: McAllen, Texas 78503  
USA

European Representative: **Unomedical a/s**  
Address: Kongevejen 2  
3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.  
The ConvaTec logo is a trademark of ConvaTec Inc.

