

# Intubation Stylets

**Intended use:**

Intubation Stylets are used as an aid to facilitate the insertion of the Tracheal Tube. The products are intended for single use only.



**Product range:**

REF*	CH size	O.D., mm (±0.2)	Length, mm (±10) (straight part)	Total length, mm (±10) (including bent part)
UM64200006	CH06	2.0	272	305
UM64200010	CH10	3.3	346	384
UM64200012	CH12	4.0	337	384
UM64200014	CH14	4.7	337	384

\* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

**Material specifications:**

**Contains latex?** ..... No  
**Contains phthalates?** ..... No  
**Stylet** .....Aluminum  
**Sheath** .....Polyvinyl chloride

**Certifications:**

**CE-marking?** ..... Yes (0123)  
**Product classification** ..... Class IIA (EU)  
**Quality system certified to** ..... ISO 13485  
**Global Medical Device Nomenclature (GMDN)** ..... 37469

**Sterilization:**

**Sterile?** .....Yes  
**Sterilized by** ..... Ethylene oxide  
**Shelf life** ..... 5 years

**Packaging specifications:**

REF	Primary package/ Single peelpack		Secondary package/ Shelf box		Tertiary package/ Shipper				Pallet
	Qty, pcs	Dimensions, mm (LxW)	Qty, pcs	Dimensions, mm (LxWxD)	Qty, pcs	Dimensions, mm (LxWxD)	Gross Wt. kg	Volume m <sup>3</sup>	Qty, pcs
UM64200006	1	500 x 66	20	510 x 80 x 45	400	520 x 340 x 250	5.2	0.044	6,400
UM64200010	1	500 x 66	20	510 x 80 x 45	400	520 x 340 x 250	6.8	0.044	6,400
UM64200012	1	500 x 66	20	510 x 80 x 45	400	520 x 340 x 250	8.0	0.044	6,400
UM64200014	1	500 x 66	20	510 x 80 x 45	400	520 x 340 x 250	9.2	0.044	6,400





## EAN codes:

REF	Primary package	Secondary package	Tertiary package
UM64200006	(01)30768455148470	(01)50768455148474	(01)70768455148478
UM64200010	(01)30768455148487	(01)50768455148481	(01)70768455148485
UM64200012	(01)30768455148685	(01)50768455148689	(01)70768455148683
UM64200014	(01)30768455148494	(01)50768455148498	(01)70768455148492

## Packaging material:

Primary package	Peelback pouch - Coated paper, Polyethylene terephthalate/Polyethylene
Secondary package	Ivory cardboard
Tertiary package	Corrugated cardboard

## Storage and handling:

	<b>FRAGILE</b> - Contents of the transport package are fragile and shall therefore be handled with care
	<b>KEEP DRY</b>
	<b>THIS WAY UP</b> - Indicates correct upright position of the transport package
	<b>DO NOT USE HAND HOOKS</b> - Hooks are prohibited for handling the transport package



**Well Lead Medical Co., Ltd.**

C-4 Jinhu Industrial Estate,  
Hualong, Panyu, Guangzhou  
511434 P.R.China

Made in China



**Shanghai Intl. Holding Corp. GmbH**

Eiffestrasse 80, 20537 Hamburg  
Germany

Distributed by



ConvaTec Limited  
First Avenue, Deeside Industrial Park  
Deeside, Flintshire, CH5 2NU, UK

[www.convatec.com](http://www.convatec.com)

**For further information on these products, please contact your local ConvaTec sales representative.**

*With reservations for misprints and alterations in product specifications. Printed version not controlled.*

*All trade marks and logos are the property of a ConvaTec group company.*

© 2017 ConvaTec Inc. AP-017959-MM, ver. JUL 2017