

Laryngeal Airway

Intended use:

The Laryngeal Airway is inserted through the mouth into the pharynx and connected to a ventilation device to secure an airway during anaesthesia. The product is intended for single use only.



List of product range:

Size	Product group	REF*	Maximum cuff inflation volume, ml	ISO connector, mm
1	Neonates/infants	MM65223010	4	15
1½	Infants	MM65223015	7	15
2	Infants/children	MM65223020	14	15
2½	Children	MM65223025	20	15
3	Children/small adults	MM65223030	30	15
4	Adults	MM65223040	40	15
5	Adults	MM65223050	50	15

List of product range:

Size	Product group	REF*	Maximum cuff inflation volume, ml	ISO connector, mm
1	Neonates/infants	MM65220010	4	15
1½	Infants	MM65220015	7	15
2	Infants/children	MM65220020	14	15
2½	Children	MM65220025	20	15
3	Children/small adults	MM65220030	30	15
4	Adults	MM65220040	40	15
5	Adults	MM65220050	50	15

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

- Contains latex?** No
- Contains DEHP?** Yes
- Tube** Polyvinyl chloride
- Inflation tube** Polyvinyl chloride
- Pilot balloon** Polyvinyl chloride
- Cuff** Polyvinyl chloride
- Connector** Polypropylene
- Biteblock** Polypropylene
- Valve** Polypropylene, Nitrile and stainless steel

Certification:

- CE-marking?** Yes (CE 0123)
- Product classification** Class IIA (EU)
- Quality system certified to** ISO 13485
- Global Medical Device Nomenclature (GMDN)**. 13364

Sterilisation:

- Sterile?** Yes
- Sterilized by** Ethylene oxide
- Shelf life** 5 years

Laryngeal Airway

Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
MM65223010	10	285 x 270 x 280	40	593 x 557 x 303	6.4	0.10	320
MM65223015	10	285 x 270 x 280	40	593 x 557 x 303	6.4	0.10	320
MM65223020	10	285 x 270 x 280	40	593 x 557 x 303	6.5	0.10	320
MM65223025	10	285 x 270 x 280	40	593 x 557 x 303	6.5	0.10	320
MM65223030	10	285 x 270 x 280	40	593 x 557 x 303	6.6	0.10	320
MM65223040	10	285 x 270 x 280	40	593 x 557 x 303	6.6	0.10	320
MM65223050	10	285 x 270 x 280	40	593 x 557 x 303	6.7	0.10	320
MM65220010	10	285 x 270 x 280	40	593 x 557 x 303	6.7	0.10	320
MM65220015	10	285 x 270 x 280	40	593 x 557 x 303	6.7	0.10	320
MM65220020	10	285 x 270 x 280	40	593 x 557 x 303	6.7	0.10	320
MM65220025	10	285 x 270 x 280	40	593 x 557 x 303	6.7	0.10	320
MM65220030	10	285 x 270 x 280	40	593 x 557 x 303	6.7	0.10	320
MM65220040	10	320 x 270 x 280	40	663 x 557 x 303	7.2	0.11	320
MM65220050	10	320 x 270 x 280	40	663 x 557 x 303	7.2	0.11	320

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Product specific information:

Biteblock

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **Unomedical Sdn Bhd**
(15290-P)
Address: Bakar Arang Industrial Estate
08000 Sungai Petani, Kedah, Malaysia

European Representative: **Unomedical a/s**
Address: Kongevejen 2
3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

