




Gas sampling lines

Intended use:

Gas sampling lines are used to deliver gas samples to a monitoring device for the measurement of respiratory gases. The product is intended for single use only.



List of product range:

Gas sampling line	REF*	Inner/outer diameter, mm	Length, m	Luer lock connectors
 895MM	895MM	1.68 / 3.33	3.0	Male/Male
	897MM	1.68 / 3.33	3.0	Male/Female
 897MM	892MM	1.27 / 3.33	3.0	Male/Male
	898MM	1.27 / 3.33	3.0	Male/Female

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
 Contains DEHP? Yes
 Gas sampling line Polyvinyl chloride
 Male/Female Acrylic

Certification:

CE-marking? Yes (CE 0120)
 Product classification Class IIA (EU)
 Quality system certified to ISO 13485
 Global Medical Device Nomenclature (GMDN) 45566

Sterilisation:

Sterile? No
 Shelf life 5 years

Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume, m ³	
895MM	*	*	30	279 x 330 x 381	1.36	2.674	750
897MM	*	*	30	279 x 330 x 381	0.91	2.674	750
892MM	*	*	30	279 x 330 x 381	0.73	2.674	750
898MM	*	*	30	279 x 330 x 381	1.36	2.674	750

* = Single packed - no shelf box

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **Unomedical Inc.**
 5701-1S, Ware Road
 Address: McAllen, Texas 78503
 USA

European Representative: **Unomedical a/s**
 Address: Kongevejen 2
 3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
 The ConvaTec logo is a trademark of ConvaTec Inc.