

Mini Vacuum Set

Wound drainage system

Intended use:

Low vacuum wound drainage systems for post operative removal of blood and body fluids. The product is intended for single use only.



Product range:

	REF*	CH/FG	Drains Shore A 69		Perforation / Introducer Length, mm
			Length, mm	I.D. / O.D., mm	
Mini vacuum set with 30 ml bellows, Redon drain with introducer and X-ray marking	1300118x	06	500	1.1 / 2.0	70 / 55
Mini vacuum set with 30 ml bellows, Ulmer drain with introducer and X-ray marking	1300218x	08	500	1.6 / 2.7	58 / 55

* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.
Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
Contains DEHP? No
Bellows Polyethylene
Drains Polyvinyl chloride
Introducer Stainless steel

Certifications:

CE-marking? Yes (CE 0086)
Product classification Class IIA (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN) 35824

Sterilization:

Sterile? Yes
Sterilized by Ethylene oxide
Shelf life 5 years

Packaging specifications:

REF	Shelf box		Shipping carton				MOQ, pcs
	Qty, pcs	Dimensions, mm (LxWxD)	Qty, pcs	Dimensions, mm (LxWxD)	Gross Wt. kg	Volume m ³	
1300118x	50	297 x 190 x 366	200	391 x 304 x 743	4.90	0.092	400
1300218x	50	297 x 190 x 366	200	391 x 304 x 743	4.90	0.092	400

Storage and handling:




ConvaTec Limited

First Avenue
 Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK

Distributed by



ConvaTec

ConvaTec Limited
 First Avenue, Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK
www.convatec.com

For further information on these products, please contact your local ConvaTec sales representative.

With reservations for misprints and alterations in product specifications. Printed version not controlled.
 All trade marks and logos are the property of a ConvaTec group company.

© 2018 ConvaTec Inc. AP-018614-MM, ver. JAN 2018