


# Siliconized Foley catheter

## Standard, 2-way, soft valve

**Intended use:**

The siliconized Foley catheter is a urinary catheter with a silicone fluid treated surface for easy passage through the urethral meatus. The siliconized catheter is intended for short-term drainage of the bladder. The product is intended for single use only.



**List of product range:**

Siliconized Foley catheter, Standard, 2-way, soft valve

CH size	5-10 ml/cc balloon	30 ml/cc balloon
12	MM51111210	MM51111230
14	MM51111410	MM51111430
16	MM51111610	MM51111630
18	MM51111810	MM51111830
20	MM51112010	MM51112030
22	MM51112210	MM51112230
24	MM51112410	MM51112430
26	MM51112610	MM51112630
28	MM51112810	MM51112830
30	MM51113010	MM51113030

**Material specifications:**

Contains latex? . . . . . Yes  
 Contains DEHP? . . . . . No  
 Catheter . . . . . Latex  
 Soft valve . . . . . Rubber  
 Sleeve . . . . . Rubber  
 Catheter surface . . . . . Silicone emulsion

**Certification:**

CE-marking? . . . . . Yes (CE 0123)  
 Product classification . . . . . Class IIB (EU)  
 Quality system certified to . . . . . ISO 13485  
 Global Medical Device Nomenclature (GMDN). . . . .34917

**Sterilisation:**

Sterile? . . . . . Yes  
 Sterilized by . . . . . Gamma radiation  
 Shelf life . . . . . 5 years

# Siliconized Foley catheter

Standard, 2-way, soft valve

## Packaging:

Siliconized Foley catheter, Standard, 2-way, soft valve

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m <sup>3</sup>	
MM51111210	10	505 x 95 x 34	800	575 x 495 x 522	23	0.16	800
MM51111230	10	506 x 95 x 34	800	576 x 495 x 522	23	0.16	800
MM51111410	10	507 x 95 x 34	800	577 x 495 x 522	23	0.16	800
MM51111430	10	508 x 95 x 34	800	578 x 495 x 522	23	0.16	800
MM51111610	10	509 x 95 x 34	800	579 x 495 x 522	23	0.16	800
MM51111630	10	510 x 95 x 34	800	580 x 495 x 522	23	0.16	800
MM51111810	10	511 x 95 x 34	800	581 x 495 x 522	23	0.16	800
MM51111830	10	512 x 95 x 34	800	582 x 495 x 522	23	0.16	800
MM51112010	10	513 x 95 x 34	800	583 x 495 x 522	23	0.16	800
MM51112030	10	514 x 95 x 34	800	584 x 495 x 522	23	0.16	800
MM51112210	10	515 x 95 x 34	800	585 x 495 x 522	23	0.16	800
MM51112230	10	516 x 95 x 34	800	586 x 495 x 522	23	0.16	800
MM51112410	10	517 x 95 x 34	800	587 x 495 x 522	23	0.16	800
MM51112430	10	518 x 95 x 34	800	588 x 495 x 522	23	0.16	800
MM51112610	10	519 x 95 x 34	800	589 x 495 x 522	23	0.16	800
MM51112630	10	520 x 95 x 34	800	590 x 495 x 522	23	0.16	800
MM51112810	10	521 x 95 x 34	800	591 x 495 x 522	23	0.16	800
MM51112830	10	522 x 95 x 34	800	592 x 495 x 522	23	0.16	800
MM51113010	10	523 x 95 x 34	800	593 x 495 x 522	23	0.16	800
MM51113030	10	524 x 95 x 34	800	594 x 495 x 522	23	0.16	800

\* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

## Storage and handling:



*With reservations for misprints and alterations in product specifications. Printed version not controlled.*

Manufacturer: **Unomedical Sdn Bhd**  
(15290-P)  
Address: Bakar Arang Industrial Estate  
08000 Sungai Petani, Kedah, Malaysia

European Representative: **Unomedical a/s**  
Address: Kongevejen 2  
3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.  
The ConvaTec logo is a trademark of ConvaTec Inc.

