




ESTEEM™  *Flex Convex*

NOW INTRODUCING OUR X-LARGE SIZE (V0) Giving you even more flexibility

The latest addition to the Esteem™+ Flex Convex range provides more options to address a wider variety of stoma and abdomen profiles

- Larger cut-to-fit range for wider diameter stomas
- up to 57mm (2 1/4")
- Helps move skin folds and creases on the abdomen out of the way
- Flexible convexity fits and moves with the contours of the body for comfort and easy positioning

 **ConvaTec** | Ostomy Care



To find out more, get in touch with one of our team members:

1-800-422-8811

FOUR CONVEXITY PLATEAU VARIANTS

The Esteem™+ Flex Convex system is available in four convexity plateau variants to accommodate a range of stomas



ESTEEM™+ FLEX CONVEX CLOSED POUCHES

Pouch size 145 x 205mm (5 11/16" x 8 1/16")

Cut-to-fit

PRODUCT CODE	PLATEAU SIZE	CUTTING RANGE	QUANTITY
421612	V1	20-43mm (4/5"-1 2/3")	30
421613	V2	20-35mm (4/5"-1 3/8")	30
421614	V3	20-25mm (4/5"-1")	30



ESTEEM™+ FLEX CONVEX UROSTOMY POUCHES

Pouch size 130 x 195mm (5 1/8" x 7 2/3")

Cut-to-fit

PRODUCT CODE	PLATEAU SIZE	CUTTING RANGE	QUANTITY
421623	V1	10-43mm (3/8"-1 2/3")	10
421626	V2	10-35mm (3/8"-1 3/8")	10
421628	V3	10-25mm (3/8"-1")	10



ESTEEM™+ FLEX CONVEX DRAINABLE POUCHES

Pouch size 150 x 220mm (5 7/8" x 8 2/3")

Cut-to-fit

PRODUCT CODE	PLATEAU SIZE	CUTTING RANGE	QUANTITY
422159	V0	20-57mm (4/5"-2 1/4")	10
421615	V1	20-43mm (4/5"-1 2/3")	10
421619	V2	20-35mm (4/5"-1 3/8")	10
421622	V3	20-25mm (4/5"-1")	10


Pre-cut

PRODUCT CODE	PLATEAU SIZE	CUTTING RANGE	QUANTITY
421624	V1	15mm (9/16")	10
421625	V1	20mm (4/5")	10
421627	V2	15mm (9/16")	10
421629	V3	15mm (9/16")	10

Pre-cut

PRODUCT CODE	PLATEAU SIZE	CUTTING RANGE	QUANTITY
421616	V1	25mm (1")	10
421617	V1	30mm (1 1/6")	10
421618	V1	35mm (1 3/8")	10
421620	V2	25mm (1")	10
421621	V2	30mm (1 1/6")	10

CONVEXITY PLATEAU VARIANTS



CONVEXITY PLATEAU WIDTH



Getting the full benefits of me+™ is easy.

Call: 1-800-422-8811 or visit www.mepluscare.com