


Oxygen Masks


Medium Concentration


Intended use:

Medium concentration masks are intended for administration of oxygen when precise control of inspired oxygen concentration is not mandatory. During inspiration the patient draws gases both from oxygen flowing into the mask through the tubing, as well as from room air via ports on the sides of the mask. Oxygen concentration delivered is 40%-60% depending on the patient's breathing. The masks are connected directly to a compressed air or oxygen supply. The products are intended for single use only.



List of product range:

Oxygen mask, Medium concentration, Elongated (under the chin) with adjustable nose-clip and elastic strap				Oxygen tubing		
	REF*	Size	Description	m	Type	Connector
	103MM	Adult	Medium concentration mask, Elongated, no tubing	-	-	-
	102MM	Adult	Medium concentration mask, Elongated, with tubing	2.1	Standard lumen	Standard
	3102MM	Adult	Medium concentration mask, Elongated, with tubing	2.1	Standard lumen	Fits-all
	9102MM	Adult	Medium concentration mask, Elongated, with tubing	2.1	No-crush lumen	Standard
	222MM	Pediatric	Medium concentration mask, Elongated, with tubing	2.1	Standard lumen	Standard
	3222MM	Pediatric	Medium concentration mask, Elongated, with tubing	2.1	Standard lumen	Fits-all
	9222MM	Pediatric	Medium concentration mask, Elongated, with tubing	2.1	No-crush lumen	Standard
	223MM	Pediatric	Medium concentration mask, Elongated, no tubing	-	-	-

Oxygen mask, Medium concentration, Standard (under the lip) with adjustable nose-clip and elastic strap				Oxygen tubing		
	REF*	Size	Description	m	Type	Connector
	203MM	Adult	Medium concentration mask, Standard, no tubing	-	-	-
	202MM	Adult	Medium concentration mask, Standard, with tubing	2.1	Standard lumen	Standard
	263MM	Infant	Medium concentration mask, Standard, no tubing	-	-	-
	3262MM	Infant	Medium concentration mask, Standard, with tubing	2.1	Standard lumen	Fits-all
	93262MM	Infant	Medium concentration mask, Standard, with tubing	2.1	No-crush lumen	Fits-all

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
 Contains DEHP? Yes
 Mask Polyvinyl chloride
 Elastic strap (latex-free) Polyester
 Nose clip Aluminum
 Mask connector/Nipple High density polyethylene
 Oxygen tubing Polyvinyl chloride
 Tube connector Polyvinyl chloride

Certification:

CE-marking? Yes (CE 0120)
 Product classification Class IIA (EU)
 Quality system certified to ISO 13485
 Global Medical Device Nomenclature (GMDN). 35171

Sterilisation:

Sterile? No
 Shelf life 5 years

Oxygen Masks

Medium Concentration

Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs	Dimensions, mm	Qty., pcs	Dimensions, mm	Gross Wt., kg	Volume, m ³	
102MM	*	*	50	406 x 355 x 254	4.26	0.037	1,250
3102MM	*	*	50	406 x 355 x 254	4.33	0.037	1,250
9102MM	*	*	50	406 x 355 x 254	4.31	0.037	1,250
103MM	*	*	50	381 x 254 x 203	2.31	0.020	1,250
202MM	*	*	50	406 x 355 x 254	3.70	0.037	1,250
203MM	*	*	50	381 x 254 x 203	1.81	0.020	1,250
222MM	*	*	50	406 x 355 x 254	3.86	0.037	1,250
3222MM	*	*	50	406 x 355 x 254	3.92	0.037	1,250
9222MM	*	*	50	406 x 355 x 254	4.08	0.037	1,250
223MM	*	*	50	381 x 254 x 203	1.81	0.020	1,250
3262MM	*	*	50	406 x 355 x 254	3.18	0.037	1,250
93262MM	*	*	50	406 x 355 x 254	4.54	0.037	1,250
263MM	*	*	50	279 x 330 x 127	0.95	0.012	1,250

* = Single packed - no shelf box

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Storage and handling:



With reservations for misprints and alterations in product specifications. Printed version not controlled.

Manufacturer: **Unomedical Inc.**
5701-1S, Ware Road
Address: McAllen, Texas 78503
USA

European Representative: **Unomedical a/s**
Address: Kongevejen 2
3460 Birkerød, Denmark

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.
The ConvaTec logo is a trademark of ConvaTec Inc.

