


neria™

Intended use:

The product is intended for the subcutaneous infusion of medication such as immunoglobulins for the treatment of Primary Immune Deficiency, deferoxamine mesylate for Thalassaemia, apomorphine for Parkinson's Disease and morphine for pain management therapy. The neria™ infusion set has been tested to be compatible with the above-mentioned drugs approved for subcutaneous infusion.


Product range:

Product code	Description	Needle gauge	Needle length	Tubing length	GTIN (Inner box)	Pcs per inner box	Pcs per outer box
neria™							
78-060-2731	neria™ infusion set	G27	10 mm	60 cm	05705244000513	10	200
78-060-2738	neria™ infusion set	G27	8 mm	60 cm	05705244000193	10	200
78-060-2736	neria™ infusion set	G27	6 mm	60 cm	05705244000476	10	200
78-080-2737	neria™ infusion set	G27	12 mm	80 cm	05705244004191	10	200
78-080-2731	neria™ infusion set	G27	10 mm	80 cm	05705244001077	10	200
78-080-2738	neria™ infusion set	G27	8 mm	80 cm	05705244000575	10	200
78-080-2736	neria™ infusion set	G27	6 mm	80 cm	05705244000551	10	200
78-110-2731	neria™ infusion set	G27	10 mm	110 cm	05705244000117	10	200
78-110-2738	neria™ infusion set	G27	8 mm	110 cm	05705244000292	10	200
78-110-2736	neria™ infusion set	G27	6 mm	110 cm	05705244000599	10	200
78-060-2931	neria™ infusion set	G29	10 mm	60 cm	05705244000216	10	200
78-060-2938	neria™ infusion set	G29	8 mm	60 cm	05705244000131	10	200
78-060-2936	neria™ infusion set	G29	6 mm	60 cm	05705244000612	10	200
78-080-2931	neria™ infusion set	G29	10 mm	80 cm	05705244000711	10	200
78-080-2938	neria™ infusion set	G29	8 mm	80 cm	05705244000674	10	200
78-080-2936	neria™ infusion set	G29	6 mm	80 cm	05705244000636	10	200
78-110-2931	neria™ infusion set	G29	10 mm	110 cm	05705244000230	10	200
78-110-2938	neria™ infusion set	G29	8 mm	110 cm	05705244000155	10	200
78-110-2936	neria™ infusion set	G29	6 mm	110 cm	05705244000759	10	200

Product specifications:

Insertion angle 90 degrees
Needle length 6 mm, 8 mm, 10 mm and 12 mm
Needle gauge 27 gauge and 29 gauge
Tubing length in cm 60, 80 and 110 cm
Tubing inner diameter G27 0,385 mm
Tubing inner diameter G29 0,305 mm
Tubing outer diameter 1,5 mm

Material specifications:

PVC free Yes
Made with natural rubber latex No
Raw materials or derivation of animal origin No
The product contains phthalates No
Introducer & connector needles may contain . . ≤ 10% nickel

Certifications:

CE-marking Yes
Product classification Class IIB (EU)
Quality system certified to ISO 13485
Tariff code HS9018390000
GMDN 35833

Sterilization:

Sterile Yes
Sterilized by Ethylene oxide
Shelf life 5 years
Usage time (drug depended) Maximum 48 hours

Drug compatibility/stability tested with¹:

Immunoglobulins
Apomorphine
Morphine (hydromorphone, morphine sulphate & morphine chloride)
Deferoxamine mesylate

Storage and handling:

	AVOID DIRECT SUNLIGHT - Keep the product away from direct sunlight and any humid environment
	KEEP DRY - Store in a dry place at room temperature
	DO NOT RE-USE
	DO NOT USE - If package is damaged
	CAUTION - See Instruction for Use



Unomedical A/S A ConvaTec Company

Aaholmvej 1-3
 Østved, 4320 Lejre, Denmark
www.neria.com